# Imuno-HAI TOXO

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## ENGLISH

### SUMMARY

Toxoplasmosis is caused by a protozoan called *Toxoplasma gondii*. As this protozoan can parasite all cells with nucleus, the infections can vary such as: coriorretinitis, myocarditis, pericarditis, pneumonia, myositis,

meningoencephalitis, hepatitis, these infections are isolated or combined, often diffuse, particularly in patients with immunodeficiency or undergoing immunosuppressive therapy.

Toxoplasmosis has emerged as a serious complication in congenitally infected children causing the death of fetus and infant or congenital disease. The symptoms such as hydrocephalus or microcephalus, coriorretinitis and cerebral calcifications are present. The transplacental congenital transmission of toxoplasmosis depends on an acute infection during pregnancy.

When the infection is acquired during the first trimester of the gestation period fetal death or major abnormalities may occur.

T. Gondii infection is insidious and characterized by fever, adynamia and adenopathy, mainly cervical. However, it is present on subclinical form presenting unpecific and light symptoms.

Toxoplasmosis diagnosis should not be made solely on the findings of one clinical assay. When making any diagnosis, it is strongly recommended to take all clinical results into consideration.

The most reliable and effective laboratorial methods for toxoplasmosis are: Indirect Immunofluorescence, passive hemagglutination and ELISA.

### PRINCIPLE OF THE METHOD

Fowl erythrocytes sensitised and stabilized with highly purified *Toxoplasma gondii* antigens show agglutination when react with antibodies against these antigens present in patient's serum or cephalo-rachidian liquor.

### KIT PRESENTATION

**REF 36096-H** (96 qualitative determinations)

- Suspension of sensitised erythrocytes with *Toxoplasma gondii* antigens (2.4ml)
- Diluent solution (1x 40ml)
- 2-Mercaptoethanol (0.5ml)
- Positive serum control (1ml)
- Negative serum control (1ml)
- Disposable "V"-well microtitration plate (1x96 wells)
- Instructions for use

**REF 36192-H** (192 qualitative determinations)

- Suspension of sensitised erythrocytes with *Toxoplasma gondii* antigens (4.8ml)
- Diluent solution (1x60ml)
- 2-Mercaptoethanol (0.5ml)
- Positive serum control (1ml)
- Negative serum control (1ml)
- Disposable "V"-well microtitration plate (2x96 wells)
- Instructions for use

**REF 36380-H** (380 qualitative determinations)

- Suspension of sensitised erythrocytes with *Toxoplasma gondii* antigens (2x4.8ml)
- Diluent solution (2x60ml)
- 2-Mercaptoethanol (1ml)
- Positive serum control (1ml)
- Negative serum control (1ml)
- Disposable "V"-well microtitration plate(4x96 wells)
- Instructions for use

### REAGENT PREPARATION AND STABILITY

• **SUSPENSION OF SENSITISED ERYTHROCYTES (1):** Ready for use. Keep at 2-8°C. Do not freeze.

Homogenize the erythrocytes prior to use by inversion and round movements, never shake the vial since it can cause haemolysis and/or auto-agglutination.

• **DILUENT SOLUTION (2):** Ready for use. Keep at 2-8°C.

• **2-MERCAPTOETANOL (3):** Protect from light. It should be diluted with the diluent solution (2). Ex: 10ml of diluent solution for 70µl of 2-mercaptoethanol (3). This mixture is intended to dilute the sample of serum in order to determine a recent infection (IgM) or the presence of non-specific antibodies.

The solution can be diluted in amber vial and stored at 2-8°C together with the kit.

• **POSITIVE SERUM CONTROL (4):** Ready for use. Keep at 2-8°C.

• **NEGATIVE SERUM CONTROL (5):** Ready for use. Keep at 2-8°C.

The kit presents good performance after being used for the first time. It is stable up to the expiration date if store at 2-8°.

### SPECIMEN COLLECTION

Fresh serum or stored at -20°C from fasting patients should be used. Freezing-thawing cycles should be avoided. Old serum can freeze in contact with 2-mercaptoethanol.

### MATERIAL REQUIRED BUT NOT PROVIDED

- Test tubes
- Serological pipettes
- Tube support and rack
- Disposable material container

### PROCEDURE

#### a. QUALITATIVE TEST ("Screening")

Objective: Intended for screening and elimination of negative serum.

- Lay the plate onto an humid piece of cloth to avoid electrostatics.
- Use 1 well of the plate for each specimen, including positive and negative controls, ready for use.

3. Dilute the serum 1:32 with diluent solution (2). Suggestion: Dilute 10µl of specimen + 310µl of diluent solution in a test tube.

4. Pipette 25µl of positive and negative controls and serum diluted at 1:32 of each sample on respective wells. Suggestion: A1 positive control, A2 negative control and A3, A4, A5... sera to be tested.

5. Add 25µl of suspension of erythrocytes (1) into each well.

6. Shake the plate by mechanical vibration (plate agitator) or beating with the fingers on the edges of the plate for 3 to 4 minutes.

7. Let rest for 1 to 2 hours at room temperature in a place free of vibrations

(IMPORTANT).

8. Read the results. (See quantitative procedure).

#### INTERPRETATION

Negative serum is regarded as NEGATIVE and mean the absence of current or recent toxoplasmosis infection. When one observes negative reaction, but one suspects toxoplasmosis, the test should be repeated after some days since pre-serological period can be present.

All sera which present positive reactions should be quantitatively tested to determine the stage of the disease

#### B. SEMIQUANTITATIVE TEST

1. Dilute the serum 1:32 according to step 3 of qualitative test.

2. Pipette 25µl of diluent solution (2) from the second well of the plate up to dilution that is intended to evaluate. Ex: If the dilution of the serum is up to 1:512, pipette 25µl of diluent in the wells A2, A3, A4 and A5. Do not pipette in the well A1.

3. Transfer 25µl of the serum diluted at 1:32 in the wells A1 and A2 of the plate. Homogenize thoroughly the serum with diluent in the well A2 (dilution 1:64) and transfer 25µl in the well A3 (dilution 1:128) and so on. Discard 25µl in the end.

4. Pipette 25µl of erythrocytes homogeneous suspension (1) in all wells.

5. Shake the plate by mechanical vibration (plate vibrator) or by beating with the fingers on the edges of the plate for 3 to 4 minutes.

6. Let rest for 1 to 2 hours at room temperature in a place free of vibration. (IMPORTANT)

7. Read the results.

#### READING

**Negative:** Erythrocytes settle in the bottom of the well forming a button.

**Positive:** Erythrocytes settle in the bottom of the well, like a rug, often presenting irregular edges.

**TITER OF THE SPECIMEN:** The highest dilution showing positive reaction. The endpoint is when the erythrocytes cover 50% of the bottom of the well.

#### DETERMINATION OF RECENT INFECTION (IgM) OR UNSPECIFIC ANTIBODIES

The use of 2-mercaptoethanol, on the hemagglutination tests, offers the possibility to differentiate the recent infection or the presence of non-specific antibodies.

#### PROCEDURE

1. Perform the test with 2-Mercaptoethanol in parallel to Quantitative Test.

2. Serum treatment: Pipette 10µl of the activated or inactivated serum + 310µl of diluted 2-Mercaptoethanol (1:32 solution) in a test tube. Incubate at 37°C for 1 hour.

#### INTERPRETATION

Negative: Patients non-immune or in pre-serological stage. In case of doubts, repeat the test after 1 or 2 weeks.

The titers of the serum with or without 2-mercaptoethanol (less than 2 dilutions): Probably recent toxoplasmosis.

**ATTENTION:** This result does not exclude recent infection (IgM) and therefore, immunofluorescence and ELISA methods for IgM detection are recommended. The drop of 2 or more titers of the treated serum related to the non-treated one: Probably recent infection (IgM)

**ATTENTION:** If titers are low and treated serum is negative, it should be avoided the presence of unpecific antibodies by repeating the test after 1 to 2 weeks and observing if the titer is higher.

**NOTE:** The hemagglutination test associated with sensitive methods such as immunofluorescence and/or ELISA enable to stabilish serological ways that best characterize the disease stages.

#### PERFORMANCE CHARACTERISTICS

In a study on 74 positive specimens performed by indirect immunofluorescence and ELISA methods, the Imuno-HAI TOXO from Wama Diagnóstica showed no false negative results which represents a sensitivity of 100%.

Another study on 268 negative specimens performed by indirect immunofluorescence and ELISA methods, the Imuno-HAI TOXO from Wama Diagnóstica showed 2 false positive results which represents a specificity of 99.3%.

#### PRECAUTIONS AND WARNINGS

1. Do not use reagents from different lots.

2. The plates should not vibrate since it can affect the reading results.

3. Erythrocytes vial should be let stand in vertical position to avoid irreversible auto-agglutination.

4. USE EXCLUSIVELY THE KIT PLATE. It should not be reused.

5. The incubation of serum with 2-mercaptoethanol in excess (more than 60 minutes) facilitates the partial freezing of serum.

6. It is acceptable variations of 1 or 2 titers for the same serum on hemagglutination reactions.

7. All human components have been tested for anti-HCV, anti-HIV and HBsAg and found negative. However, this does not assure the absence of these and other infectious diseases. All human components should be handled as potentially hazardous.

8. Do not freeze reagents.

9. Sodium azide 0.1% is used as preservative. It can react with lead and copper plumbing to form highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide buildup. Sodium azide is toxic if ingested.

10. Disposal in accordance with local regulations.

11. Follow good laboratory practices (GLP) related to storage, handling and material disposal.

#### WARRANTY

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be

invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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