MS 10310030069
10310030090*

Imuno-Látex ASLO

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ENGLISH

SUMMARY

Streptolysin-O, produced by almost all strains of *Streptococcus pyogenes*, is an haemolysin oxygen-labile (inactivated by oxygen) and it is highly antigenic.

Streptolysin-O has antigenic activity which results in the production of antibodies known as anti-streptolysin-O. These antibodies have been universally accepted as diagnostic procedure in streptococcal infections such as rheumatic fever and acute glomerulonephritis. Anti-streptolysin-O antibodies concentration increases at the end of the first week after streptococcal infection and reaches its highest value between third and fifth weeks. The concentration has its normal values from second to fourth month.

PRINCIPLE OF THE METHOD

Latex particles purified and stabilised coated with streptolysin-O show clear agglutination when mixed, in a test card, with serum containing high levels of anti-streptolysin antibodies.

KIT PRESENTATION

REF 21060-L (60 determinations)

1. Suspension of latex coated with streptolysin-O (1,5ml)
2. Positive control serum (0,5ml)
3. Negative control serum (1ml)
4. Stirrers (60)
5. Test cards (2)
6. Instructions for use

REF 210100-L (100 determinations)

1. Suspension of latex coated with streptolysin-O (2,5ml)
2. Positive control serum (0,5ml)
3. Negative control serum (1ml)
4. Stirrers (100)
5. Test cards (2)
6. Instructions for use

REF 21000-L (2.5ml Latex)

1. Suspension of latex coated with streptolysin-O (2.5ml)
2. Instructions for use

MATERIAL REQUIRED BUT NOT PROVIDED

- Test tubes
- Test tube shelf rack for tips
- Serological pipettes
- Disposable material container
- Saline 0,9%

REAGENT STABILITY AND STORAGE

• **SUSPENSION OF LATEX (1):** Stable if stored at 2-8°C up to expiration date. Allow the latex to reach room temperature and mix it prior to use. Do not freeze. It contains Sodium azide 0.1%.

• **POSITIVE CONTROL SERUM (2):** Ready for use. Stable if stored at 2-8°C up to expiration date. Allow the control to reach room temperature prior to use. Anti-streptolysin-O concentration: > or = 200UI/ml. It contains Sodium azide 0.1%.

• **NEGATIVE SERUM CONTROL (3):** Ready for use. Stable if stored at 2-8°C up to expiration date. Allow the control to reach room temperature prior to use. It contains Sodium azide 0.1%.

The kit presents good performance after being used for the first time. It is stable up to the expiration date if store at 2 - 8°C.

SPECIMEN COLLECTION AND STORAGE

Hemolytic, lipaemic or contaminated serum must be discarded. Specimen may be stored at 2-8°C for up to 48 hours prior testing. If longer storage is required, store in freezer at - 20°C for up to 6 weeks. Do not use diluted serum. Do not use plasma since fibrinogen can form non-specific agglutination.

PROCEDURE

A. Qualitative test

1. Pipette 25µl of patient's serum in the area of test card.
2. Gently shake the suspension of latex (1) and pipette 25µl in the same area of test card.
3. Using the stirrers, mix the serum dilution and the latex and spread them over the entire circle.
4. Rotate the test card for 2 minutes and observe the agglutination while holding the slide under a good source of light.

ATTENTION: Positive and negative controls should be included at regular intervals in order to check the procedure and reagents.

Prozone effect can be present in concentrations higher than 1500IU/ml. If it is the case, the specimen should be diluted.

READING

Positive:

Clear agglutination is seen. Concentration equal or > 200IU/ml.

Negative: No agglutination is seen. Concentration <200IU/ml.

ATTENTION: The sensitivity was adjusted to detect 200IU/ml. Therefore, negative serum is lower than 200IU/ml.

B. Semi-Quantitative test

1. Using saline (NaCl 0.9%) prepare serial dilutions of the patient's serum at 1:2, 1:4, 1:8, 1:16 and so on.
2. Pipette 25µl of each dilution in each area of test card.
3. Mix the suspension of latex (1) and pipette 25µl in each area where the sample dilution is present.
4. Using the stirrer, mix the serum dilution with latex and spread over the entire circle (Use a separate stirrer for each dilution).
5. Rotate the test card for 2 minutes and observe the agglutination while holding the slide under a good source of light. (See Reading).

The highest serum dilution showing a definite agglutination pattern is considered the titer end point. ASO concentration is calculated as follow: Concentration (IU/ml) = 200 x D, where 200 is the sensitivity of test and D is the highest dilution showing agglutination.

EXAMPLE: If the titer is 1:8, concentration of ASO on the specimen is: 200x8 = 1600IU/ml.

INTERPRETATION

ASO normal values can vary with the age of the person. Normally this value is below 250IU/ml. Specialists consider 333IU/ml as normal values up to 5 years of old and up to 500 IU/ml values above 5 years of old.

It is important to highlight that ASO does not produce much information about the current state of the disease. The levels of the titers and severity of the symptoms are not necessarily parallel.

ASO levels can remain high with fixed values for months even with normal degrees of infection.

Early intake of antibiotics and corticosteroid up to 12 to 15 days after streptococcal infection may decrease the immunological response and consequently show ASO with low titers.

Although ASO is regarded as an important test, it is only a previous response to streptococcal infection and it is not a cause of rheumatism.

PRECAUTIONS AND WARNINGS

1. Store the reagents at 2-8°C. Do not freeze.
2. Read the results after 2 minutes can cause false positive results.
3. At the end of the test, wash the test card with distilled water. If prefer to wash later, use water with a neutral detergent and rinse out slightly with distilled or deionized water. Dry prior to use. Vestigium of detergent may affect the results by causing false positive results.
4. Use separate stirrer for each sample.

5. **Imuno Látex ASLO** (latex and controls) contain 0.1% sodium azide as a preservative which may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. On disposal, flush with large quantities of water.

6. The source material used in the preparation of the controls has been tested and confirmed negative for HCV, HIV and HBSAg. Because no test can offer complete assurance that products derived from human source will not transmit infectious agents it is recommended that these controls be treated as potential biohazards materials.

7. Disposal in accordance with local regulations.

8. Follow the good laboratory practices (GLP) related to storage, handling and material disposal.

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