Imuno-Rápido HCV

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ENGLISH

SUMMARY

Hepatitis C Virus (HCV) is known to be the major cause of parenterally transmitted non-A, non-B hepatitis. HCV is a small, enveloped, single-stranded RNA virus of the Flaviviridae family. The incubation period is from 6 to 7 weeks, with an average of 2 to 26 weeks. Approximately 20 to 30% of patients with acute infection present jaundice, 70% of patients are asymptomatic, and 70 to 85% of them develop chronic hepatitis as well as cirrhosis and carcinoma hepatocellular.

The risk factors associated with the transmission of HCV are percutaneous exposure to blood through blood transfusion and infected blood donors, intravenous drug use and kidney dialysis and accidents with cuts or perforations.

Other risks include sexual intercourse with infected partners and contact with positive anti-HCV people. Associated to these risk factors, it is found a high percentage (45%) of patients with no specific factor who belong to a low socio-economic level.

The development of diagnostic tests for non-A and non-B hepatitis is linked to cDNA

cloning obtained from HCV genome by Qui-Lim Choo et al., in 1989. This genome has around 9,400 nucleotides and at least 9 codified proteins, 3 structural (Core, E1 and E2) and 6 non-structural (NS2, NS3, NS4A, NS4B, NS5A and NS5B). The test device from Imuno-Rápido HCV has a mixture of synthetic and recombinant proteins from HCV.

PRINCIPLE OF METHOD

Anti-HCV antibodies present in the specimens bind to the colloidal gold conjugate and migrate upward on the membrane of the test device by capillary action to react with the recombinant HCV antigen immobilized in the positive reaction area (T) generating a pink colored band. The absence of anti-HCV antibodies indicates a negative result. The remaining complex continues to migrate to a control area in the test device and produces a pink colored band in the control area (C). This colored band indicates that the test has been performed properly.

KIT PRESENTATION

REF 621010-R – 10 determinations

1. Test device: 10 units

2. Diluent solution: 1x2 ml

3. Instructions for use.

REF 621020-R – 20 determinations

1. Test device: 20 units

2. Diluent solution: 2x2 ml

3. Instructions for use.

REF 621040-R – 40 determinations

1. Test device: 40 units

2. Diluent solution: 4x2 ml

3. Instructions for use.

MATERIAL REQUIRED BUT NOT PROVIDED

Paper towel

Disposable material container

Automatic pipette

REAGENT STABILITY AND STORAGE

TEST DEVICE (1): Store at 2-30°C. Do not freeze.

DILUENT SOLUTION (2): Ready for use. Stable if stored at 2-30°C up to expiration date. It contains Sodium azide 0.1%. Do not freeze.

Allow the test device, specimen and diluent solution to equilibrate to room temperature prior to testing.

SPECIMEN COLLECTION AND HANDLING

This test can be run on serum or whole blood.

Fresh serum specimen is preferable. Do not use haemolysed, clotted, contaminated, viscous/turbidity and lipaemic specimen. Testing should be performed immediately after specimens have been collected. Do not leave the specimen at room temperature for prolonged periods. The specimen may be stored at 2-8°C up to 48 hours. For longer term storage, serum specimens should be kept below -20°C.


Whole blood specimen should be collected by fingerstick without anticoagulants since the presence of low reactivity specimens may cause false negative results. Specimens should not be frozen and thawed repeatedly.

ATTENTION: Frozen specimens must be completely thawed and mixed well. Keep the specimens in a vertical position for complete particles sedimentation. Do not shake the specimen.

Diluted specimens can yield false negative results.

PROCEDURE

1. Allow test device (1) to equilibrate to room temperature prior to removing it from the pouch.
2. Dispense 10 µl of specimen (avoid trapping air bubbles) in the specimen well of the test device.


1. Dispense 3 drops (100 µl) of diluent solution (2).


2. Test should be read from 10 to 15 minutes. Do not interpret the results after 20 minutes.

INTERPRETATION OF RESULTS

NEGATIVE: One pink colored band appears in the control area (C).



POSITIVE: Two pink colored bands appear. One in the test area (T) and one in the control area (C).

Obs. Interpret as positive any intensity of the pink band in the test area.

INVALID: No band appears in the test (T) and control (C) areas or no band appears in the control area (C).

Obs: Do not interpret the results after the prescribed time limit.

PERFORMANCE CHARACTERISTICS

Immuno-Rápido HCV from WAMA Diagnóstica showed 100% of sensitivity using 127 known positive specimens. Among 777 known negative specimens, it was found 01 false positive result. The results show that the sensitivity of the Immuno-Rápido HCV from Wama Diagnóstica is 99.8%.

LIMITATIONS

Immuno-Rápido HCV is used for the detection of anti-HCV antibodies in serum and whole blood specimens, therefore, repeatedly positive results with this test is an evidence of antibodies in the specimens. Positive results should be confirmed by RIBA (Recombinant Immunoblot Assay) or another confirmatory test.

False positive and false negative results may be yielded with this test whose proportion will depend on the prevalence of the disease on the assayed population.

As with all diagnostic tests, all results must be considered with other clinical information available to the physician.

PRECAUTIONS AND WARNINGS

1. For professional "in vitro" diagnostic use only.
2. Read carefully this instructions for use insert prior to testing.
3. The expiration date is the last day of the month printed on test device pouch and on the kit box.
4. Keep the kit away from heat and sunlight.
5. Do not freeze the test device. Irreversible deterioration may be yielded.
6. Sodium azide 0.1% may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide might be toxic if ingested.
7. Allow the kit components to equilibrate to room temperature prior to testing.
8. Do not use the kit beyond the expiration date.
9. Do not interchange kit components with those from sources other than the same catalog number from WAMA DIAGNOSTICA.
10. Disposal in accordance with local regulations.
11. Follow the good laboratory practices (GLP) in storing, dispensing and disposing of these materials.

WARRANTY

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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