MS 10310030070  
10310030097\*

# Waaleroose

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## ENGLISH

## SUMMARY

Rheumatoid arthritis diagnosis is based on clinical assay but the

laboratorial and radiological tests are useful to confirm the clinical diagnosis and to evaluate the seriousness and the phase of the disease.

Rheumatoid Factor (RF) in serum is one of the main clinical markers for rheumatoid arthritis. RF is a term used to describe a group of antibodies (IgM,IgG,IgA and IgE) directed against the antigenic site on the Fc portion of the IgG antibodies. Many tests have been suggested for the detection of RF such as tests which use human or sheep erythrocytes coated with human or animal IgG. Other methods have also been described such as latex particles coated with human IgG. **Waaleroose** is a very sensitive test using sheep red blood cells coated with purified and stabilized rabbit IgG.

## PRINCIPLE OF THE METHOD

When the suspension of erythrocytes is mixed in the circle of the test card with serum containing high levels of rheumatoid factor, a clear agglutination can be seen in 2 minutes.

## KIT PRESENTATION

**REF 5860-A (60 determinations)**

1. Suspension of sheep erythrocytes coated with rabbit IgG anti-sheep erythrocytes (1x1.5ml)
2. Positive serum control (1x0.5ml)
3. Negative serum control (1x1.0ml)
4. Stirrers (60)
5. Test card (2)
6. Package insert

**REF 58100-A (100 determinations)**

1. Suspension of sheep erythrocytes coated with rabbit IgG anti-sheep erythrocytes (1x2.5ml)
2. Positive serum control (1x0.5ml)
3. Negative serum control (1x1.0ml)
4. Stirrers (100)
5. Test card (2)
6. Package insert

**REF 5800-A (2.5 ml)**

1. Suspension of sheep erythrocytes coated with rabbit IgG anti-sheep erythrocytes (1x2.5ml)

2. Package insert.

## REAGENT STABILITY AND STORAGE

**SUSPENSION OF SHEEP ERYTHROCYTES (1):** Bring the suspension to room temperature and mix gently prior to use. Stable if stored at 2-8°C up to expiration date. **Do not freeze.** It contains sodium azide 0.095%.

**POSITIVE CONTROL SERUM (6):** Ready for use. Stable if stored at 2-8°C up to expiration date. It contains Sodium azide 0.095%.

**NEGATIVE SERUM CONTROL:** Ready for use. Control should be brought to room temperature prior to use. Stable if stored at 2-8°C up to expiration date. It contains sodium azide 0.095%. **Note.** The kit retains the same performance after first use and is stable until the expiration date described on the label, provided it is kept at the indicated temperature (2-8 °C).

## SPECIMEN COLLECTION AND HANDLING

Only serum specimen should be used with this test. Do not use haemolysed, lipaemic or contaminated serum. Specimen should be stored at 2-8°C for up to 48 hours prior to test. If longer storage is required, store at - 20°C for up to 6 weeks. Undiluted serum specimen should be used. Do not use plasma as this will adversely affect the results (non-specific agglutination).

## MATERIAL REQUIRED BUT NOT PROVIDED

- Test tubes
- Serological pipettes
- Test tube rack and tip rack
- Disposable material container
- Saline 0.9%

## PROCEDURE

### A. Qualitative test

1. Bring the reagents and patient's serum to room temperature prior to use.
2. Pipette 25µl of patient's serum in the circle of the test card.
3. Gently shake the suspension of sheep erythrocytes (1) and pipette 25µl in the same circle of the test card.
4. Using the stirrers, mix the serum with the suspension and spread

them over the entire circle.  
5. Leave the test card on a flat surface for 2 minutes.  
6. Immediately after that time, tilt the test card once to about 45° and again leave the test card on a flat surface for 1 minute.  
7. Observe the test card for any degree of agglutination by GENTLY tilting the test card (**important**).  
**ATTENTION:** Positive and negative controls should be run daily in order to check the optimal reactivity of the reagent.

#### READING

**Positive:** Indicated by a clear degree of agglutination. Concentration equal or > 8IU/ml.  
**Negative:** Indicated by no visible agglutination pattern. Concentration <8IU/ml.

**ATTENTION:** Sensitivity is 8IU/ml according to WHO international standard.

#### B. Semi-Quantitative test

1. Dilute patient's serum with saline (Na Cl 0.9%) 1:2, 1:4, 1:8, 1:16 and so on.  
2. Pipette 25µl of each dilution in each circle of the test card.  
3. Gently shake the suspension of sheep erythrocytes (1) and pipette 25µl in the same circle where the specimen diluted is present.  
4. Using the stirrers, mix the serum with the suspension and spread them over the entire circle. (Use a separate stirrer for each dilution).  
5. Place the test card on a flat surface for 2 minutes.  
6. Immediately after that time, tilt the test card once to about 45° and again leave the test card on a flat surface for 1 minute.  
7. Observe the test card for any degree of agglutination by GENTLY tilting the test card (important).

The highest dilution of the serum which shows a clear degree of agglutination is considered the titer end point. The concentration of RF is calculated as follow:  
Concentration (IU/ml) = 8 x D, where 8 is the sensitivity of the test and D is the highest dilution which shows agglutination.

**EXAMPLE:** If the titer is 1:8, the concentration of RF on specimen is: 8x8 = 64IU/ml.

#### INTERPRETATION

Of those patients with clinical diagnosis of rheumatoid arthritis approximately 70-80% are seropositive for rheumatoid facto. RF is not specific of arthritis because patients suffering from syphilis, LES, infectious mononucleosis, hepatitis, hypergammaglobulinaemia etc, can have positive results diagnosis, but with low titers.

A positive result can be expected in less than 5% of healthy individuals.

RF positive result is present 6 months from the onset of the disease.

#### TEST PERFORMANCE

Waller Rose Kit presented 100% of sensibility and specificity. It was used 50 specimens of known results and they were compared against other commercial kit.

#### PRECAUTIONS AND WARNINGS

1. Store the reagents at 2-8°C. **Do not freeze.**  
2. The sensitivity of the test may be reduced at low temperatures. The best results are achieved above 10°C.  
3. Delays in reading the results may generate false positive results.  
4. The results obtained with **Waller Rose** do not compare with those of **Imuno-Latex FR**. The difference of the results between these methods do not reflect the differences of both to detect rheumatoid factors.  
5. At the end of the test, wash the test card with distilled water. If it is not possible, use water with a neutral detergent and rinse out exhaustively with distilled or deionized water. Dry the test card prior to use.  
6. Use separate stirrer for each sample.  
7. Waaler Rose ( red blood cells and controls) contain 0.095% sodium azide as a preservative which may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. On disposal, flush with large quantities of water.  
8. The source material used in the preparation of the controls has been tested and confirmed negative for HCV, HIV and HBsAg.

Because no test can offer complete assurance that products derived from human source will not transmit infectious agents, it is recommended that these controls be treated as potential biohazards materials.

#### WARRANT

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replace if technical support finds evidence that running, handling and storage were not properly followed.

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