Imuno-CON TOXO

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ENGLISH

SUMMARY
Toxoplasmosis is caused by a protozoan called *Toxoplasma gondii*. As this protozoan can parasite all cells with nucleus, the infections can vary such as: coriorretinitis, myocarditis, pericarditis, pneumonia, myositis,

meningoencephalitis, hepatitis. These infections are isolated or combined, often diffuse, they can be seen in patients with immunodeficiency or undergoing immunosuppressive therapy.

Toxoplasmosis has emerged as a serious complication in congenitally infected children causing the death of fetus and infant or congenital disease. The symptoms such as hydrocephalus or microcephalus, coriorretinitis and cerebral calcifications are present. The transplacental congenital transmission of toxoplasmosis depends on an acute infection during pregnancy

When the pregnancy is acquired during the first trimester of the gestational fetal, death or major abnormalities may occur.

T. gondii infection is insidious and characterized by fever, adynamia and adenopathy, mainly cervical. However, it is present on subclinical form presenting non-specific and light symptoms.

Toxoplasmosis diagnosis should not be made solely on the findings of one clinical assay. When making any diagnosis, it is strongly recommended to take all clinical results into consideration.

The most reliable and effective laboratorial methods for toxoplasmosis are: Indirect. Immunofluorescence, passive hemagglutination and ELISA.

PRINCIPLE OF THE METHOD

Anti-toxoplasmosis antibodies present in the serum bind to an antigen-fixed slide and are revealed by an antigammaglobuin labelled with FITC.

KIT PRESENTATION

[REF 1660-I (60 determinations)

- 6-well *Toxoplasma gondii* Slides (10 slides)
- Antigammaglobulin G human labelled with FITC (1x4ml)
- Phosphate Buffered Saline (PBS) 20x concentrated (2x50ml).
- Mounting Medium (4 ml)
- Evans Blue (2 ml)
- Human positive control serum (1 ml)
- Human negative control serum (1 ml)
- Coverslips (10 units)
- Instructions for use

[REF 16100-I (100 determinations)

- 10-well *Toxoplasma gondii* Slides (10 slides)
- Antigammaglobulin G human labelled with FITC (1x5ml)
- Phosphate Buffered Saline (PBS) 20x concentrated (3x50ml).
- Mounting Medium (4 ml)
- Evans Blue (2 ml)
- Human positive control serum (1 ml)
- Human negative control serum (1 ml)
- Coverslips (10 unidades)
- Instructions for use

REAGENT STABILITY AND STORAGE

• **SLIDES (1):** Let the pouch equilibrate to room temperature for 15 minutes then remove the slide from the pouch. Stable if stored in refrigerator at 2°C to 8°C up to expiration date.

• **CONJUGATE (2):** Ready for use. Stable if stored at 2-8°C up to expiration date. Protect from light. It contains Sodium Azide 0.1%.

• **PHOSPHATE BUFFERED SALINE (PBS) (3):** Dilute PBS at 1:20 (Ex: 10ml of PBS + 190ml of distilled or deionized water). Stable up to expiration date if stored at 2-8°C in a clean and sealed container. Discard if turbidity or changing on pH is present.

Obs.: In case of PBS crystallization is present before dilution, keep it in water bath at 37°C up to a complete crystal dissolution.

• **MOUNTING MEDIUM (4):** Ready for use. Stable if stored at 2-8°C up to expiration date.

• **EVANS BLUE(5):** Ready for use. Stable if stored at 2-8°C up to expiration date.

• **POSITIVE CONTROL SERUM (6):** Ready for use. Stable if stored at 2-8°C up to expiration date. It contains Sodium azide 0.1%.

• **NEGATIVE CONTROL SERUM (7):** Ready for use. Stable if stored at 2-8°C up to expiration date. It contains Sodium azide 0.1%.

The kit presents good performance after begin used for the first time. It is stable up to the expiration date if store at 2-8°C.

SPECIMEN COLLECTION AND STORAGE

Only serum specimens should be used for this procedure. Grossly haemolysed, lipaemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store serum at 2-8°C up to 72 hours. For longer storage, the serum should be stored in freezer at -20°C. Avoid repeated thawing and freezing of specimens.

MATERIAL REQUIRED BUT NOT PROVIDED

- Fluorescence microscope
- Micropipette
- Test tube
- Serological pipettes
- Coplin jar or similar
- Distilled or deionized water
- 1 liter container
- Incubation chamber
- Paper towel

PROCEDURE

QUALITATIVE TEST

Objective: It's intended to screening and removal of negative serum.

- Dilute an unknown serum at 1:32. Suggestion: 0.31ml of PBS with 10µl of specimen.

- Use undiluted control serum.
- Follow the protocol from step 2 of semi-quantitative test.

SEMI-QUANTITATIVE TEST

Objective: Positive serum in the qualitative test are titrated in order to determine the stage of the disease.

- Dilute positive specimens with PBS, starting from a dilution at 1:32, following the ratio 2.

Ex: 1:32, 1:64, 1:128, 1:256..., 1:4000 and so on.

- Allow the pouch to equilibrate room temperature for 15 minutes, then remove the slide from the pouch. Remove carefully the slide without touching on the substrate.

- Add 1 drop of positive control (6) and negative (7), undiluted, in the wells 1 and 2 respectively. Avoid overfilling the wells.

- Add 1 drop of unknown serum diluted at 1:32 in the qualitative test and 1 drop of titrated serum in the semi-quantitative test in the remaining wells. Avoid overfilling the wells.

- Incubate the slide for 30 minutes at room temperature.

- Remove the slides from the incubation chamber. Hold the slide at tab end and rinse it with 10ml (PBS) (3). With a pipette, lead PBS by the longitudinal edge of the slide. Be sure not to touch reagent wells to avoid damaging the substrate. Transfer the slide into the Coplin jar or similar and rinse it three times with PBS for 5 minutes each time. Gently shake the jar during rinsing.

- Remove the slide from Coplin jar, one by one. Blot the edge of slide on a paper towel to remove PBS. Dry around reactive wells. To prevent the slide from drying, go IMMEDIATELY to step 8.

- Place the slide back to incubation chamber. Apply 1 drop of the conjugate in each well. Be sure to fill the well completely.

- Incubate the slide for 30 minutes at room temperature. Protect it from light.

- Remove the slide from incubation chamber. Hold the slide at tab end and rinse it with 10 ml of PBS(3). With a pipette, lead PBS by the longitudinal edge of the slide. Be sure not to touch reagent wells to avoid damaging the substrate. Transfer the slide into Coplin jar or similar and rinse it three times with PBS for 5 minutes each time. Gently shake the jar. Add 2 to 3 drops of Evans Blue (5) in the Coplin in the last rinsing.

- Remove the slide from Coplin jar, one by one. Blot the edge of the slide on a paper towel to remove PBS excess. Dry around reactive wells. To prevent wells from drying, go IMMEDIATELY to step 12

- Apply 3 to 4 drops of Mounting Medium (4) in the reactive wells. Cover the slides with coverslips avoiding air bubbles. Dry the excess of Mounting Medium with a paper towel. Clean the slide reverse.

- The slide must be read under a fluorescence microscope as soon as possible. If it is not possible, store the slide at 2-8°C, protect it from light and interpret in the following day. Mounting Medium should not dry, if necessary, add more.

READING

Negative: ABSENCE of yellow-greenish fluorescence on all toxoplasma outline. Parasites are seen by a redish color.

Positive: PRESENCE of a typical yellow-greenish fluorescence on all toxoplasma outline (antibody "anti-wall").

OBSERVATION: Always check positive and negative controls for reaction control purposes.

INTERPRETATION

The interpretation of indirect immunofluorescence results to evaluate the absence of a current or past infection should be based on the presence or absence of anti-toxoplasma antibodies IgG and IgM.

Serum negative for IgG and IgM anti-toxoplasma are regarded as NEGATIVE and mean the absence of current or past toxoplasmosis infection.

Serum positive for IgG (mostly high titers) and IgM anti-toxoplasma are regarded as POSITIVE (Ex: IgG = POSITIVE up to 1:8000 and IgM = POSITIVE up to up to 1:256) and mean recent toxoplasma infection.

Serum positive for IgG and negative for IgM anti-toxoplasma mean past infection.

PRECAUTIONS AND WARNINGS

- For in vitro Diagnostic use only.
- The kit should be stored at 2-8°C.
- The conservation of the slide is done by sealing the coverslip edge with nail varnish and storing at 2-8°C in a dark room for 2 months.

- Sodium Azide 0.1% is used and it can react with lead and cooper plumbing forming highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide buildup. Sodium azide is toxic if ingested.

- All human components have been tested for HBsAg, anti-HIV, anti-HCV and found to be negative. However, this does not assure the absence of these and other infectious diseases. All human components should be handled as potentially hazardous.

- Do not interchange kit components with those from sources other than the same code number from WAMADIAGNOSTICA.

- Do not use components beyond their expiration date.

- Disposal in accordance with local regulations.

- The periodical maintenance of the microscope is recommended since an expired lamp can interfere in the result analysis.

- Follow good laboratory practices (GLP) related to storage, dispensing and material disposal.

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