MS 10310030077

Imuno-CON CHAGAS

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ENGLISH

SUMMARY

Chagas disease or American Trypanosomiasis is an endemic and chronic disease. Its pathogenic agent is a flagellate protozoan called *Trypanosoma cruzi* which is transmitted to humans by a bug, triatominae.

3 stages can be present: acute, latent or indetermined and chronic.

The acute phase can present manifestations that include fever, myocardiopathy, lymphadenopathy, hepatoesplenomegaly and parasitemia. Pseudocyst is formed after the intracellular multiplication of parasite on muscles and on cells of endothelial reticulum system.

Latent and intermediary stages are asymptomatic. The disease can evolve for

the chronic phase with manifestations of myocardiopathy, degeneration of ganglion cells of peripheral and central nervous system, hypertrophy and dilatation of some organs such as oesophagus and colon.

Chagas disease remains a major health problem in many Latin American countries because the high numbers of deaths.

As the minority of patients presenting positive sorology for *T. cruzi* develop clinical evidence for chronic disease, information provided by the laboratory has a decisive and fundamental role for etiologic diagnostic.

Chagas disease methods are: Complement fixation, agglutination, precipitation, immunofluorescence, hemagglutination and ELISA from which the most used are indirect hemagglutination (HA), indirect immunofluorescence (IFI) and ELISA.

PRINCIPLE OF THE METHOD

Antibodies anti-*Trypanosoma cruzi* present in the serum bind to an antigen fixed on slide, a fluorescein-labeled, anti-human IgG conjugate (anti-human globulin labelled with FITC) is used to visualize the reaction.

KIT PRESENTATION

REF 1460-I (60 determinations)

- 6-well *Trypanosoma cruzi* slide (10 slides)
- Fluorescein-labelled, anti-human IgG conjugate (5ml)
- Phosphate Buffered Saline (PBS) 20x (3 x 50ml)
- Mounting Medium (4ml)
- Evans Blue (2ml)
- Positive control serum (1ml)
- Negative control serum (1ml)
- Coverslips (10 units)
- Instructions for use

REF 14100-I (100 determinations)

- 6-well *Trypanosoma cruzi* slide (10 slides)
- Fluorescein-labelled, anti-human IgG conjugate (5ml)
- Phosphate Buffered Saline (PBS) 20x (3 x 50ml)
- Mounting Medium (4ml)
- Evans Blue (2ml)
- Positive control serum (1ml)
- Negative control serum (1ml)
- Coverslips (10 units)
- Instructions for use

REF 14200-I (200 determinations)

- 10-well *Trypanosoma cruzi* slides (20 slides)
- Fluorescein-labelled, anti-human IgG conjugate (2x5ml)
- Phosphate Buffered Saline (PBS) 20x (6 x 50ml)
- Mounting Medium (2x4ml)
- Evans Blue (2x2ml)
- Positive serum control (2x1ml)
- Negative serum control (2x1ml)
- Coverslips (20 units)
- Instructions for use

STABILITY AND STORAGE

• TRYPANOSOMA CRUZI SLIDE (1):

Let pouch equilibrate room temperature, then remove slide from the pouch. Stable if stored at 2-8°C up to expiration date.

• **FLUORESCIN-LABELED, ANTI-HUMAN IgG CONJUGATE (2):** Ready for use. Stable if stored at 2-8°C up to expiration date. Contains sodium azide 0,095%.

• BUFFERED SALINE (PBS) (3):

Dilute PBS 1:20 (Ex.: 10ml of PBS + 190ml of distilled or deionized water). Stable up to expiration date if stored at 2-8°C in a clean and sealed container. Discard if turbidity or changing on pH is present.

Obs.: In case PBS crystallization is present before the dilution, keep it in water bath at 37°C up to a complete crystal dissolution.

• **MOUNTING MEDIUM (4):** Ready for use. Stable if stored at 2-8°C up to expiration date. Contains sodium azide 0,095%.

• **EVANS BLUE (5):** Ready for use. Stable if stored at 2-8°C up to expiration date.

• **POSITIVE CONTROL SERUM (6):** Ready for use. Stable up to expiration date.

• **NEGATIVE SERUM CONTROL (7):** Ready for use. Stable up to expiration date.

Obs.: Kit presents results after its first handling and it is stable up to expiration date if stored at 2-8°C.

SPECIMEN COLLECTION

Only serum specimens should be used for this procedure. Grossly haemolysed, lipemic or microbially contaminated specimens may interfere in the performance of this test and should not be used. Store serum at 2-8°C up to 72 hours. For long time storage, serum should be stored at -20°C. Avoid repeated thawing and freezing cycles.

MATERIAL REQUIRED BUT NOT PROVIDED

- Fluorescence microscope
- Serological pipette
- Coplin jar or similar
- Test tubes and rack
- Distilled or deionized water
- 1 liter container
- Incubation chamber
- Paper towel

PROCEDURE

- Inactivate serum for 30 minutes at 56°C.

Unknown serum 1:30 (suggestion: 50µl of serum + 1.45ml of PBS).
Pour the pouch equilibrated to room temperature for 15 minutes, then remove the slide from the pouch. Carefully remove the slide without touching on the slide. Label it and store into moisture chamber.

1 drop of positive control(6) and negative (7), **undiluted**, in the control Avoid overfilling the wells.

1 drop (50µl) of unknown serum, 1:30 dilution, in the remaining wells. overfilling the wells.

Substrate slide for 30 minutes at room temperature. Hold the slide at tab end and with 10ml (PBS) (3). With a pipette, lead PBS by the longitudinal edge of the slide. Be sure not to touch on the reagent wells to avoid damaging the slide. Transfer the slide into Coplin jar or similar and rinse it three times with PBS for 5 minutes each time. Gently shake the jar during rinsing.

Remove the slide from Coplin jar, one by one. Blot the edge of the slide on a towel to remove PBS. Dry around reactive wells. To prevent slide from going IMMEDIATELY to step 9.

Remove the slide back to incubation chamber. Apply 1 drop of Fluorescein-anti-human IgG conjugate on each well. Be sure to fill the well completely.

Substrate slide for 30 minutes at room temperature. Protect from light. Remove the slide from incubation chamber. Hold the slide at tab end and with 10 ml of PBS(3). With a pipette, lead PBS by the longitudinal edge of the slide. Be sure not to touch on the reagent wells to avoid damaging the slide. Transfer the slide into Coplin jar or similar and rinse it three times with PBS for 5 minutes each time. Gently shake the jar. Add 2 to 3 drops of Evans Blue in Coplin on the last rising.

Remove slide from Coplin jar, one by one. Blot the edge of the slide on a towel to remove PBS excess. Dry around reactive wells. To prevent from drying, go IMMEDIATELY to step 13.

Apply 3 to 4 drops of Mounting Medium (4) on reactive wells. Cover slips with tape avoiding bubbles. Dry Mounting Medium excess with a paper towel. The slide reverse.

The slide must be read under a fluorescence microscope as soon as possible. However, in case the reading is not possible, store the slide at 2-8°C, protected from light and read the following day. Mounting Medium should not dry, otherwise, apply more.

READING
Positive: Presence of yellow-greenish fluorescence on all Trypanosoma cruzi outline. Parasites are seen by a redish color.

NEGATIVE
Absence of a typical yellow-greenish fluorescence on all Trypanosoma cruzi outline.

CAUTION: Always check positive and negative controls for reaction purposes.

PRETESTING
A positive/negative reaction is regarded as NEGATIVE, meaning the patient has no parasites. Serum presenting positive reaction is regarded as POSITIVE, meaning patients with Trypanosoma cruzi parasites.

No diagnosis should be based on a single serologic test, including serologic tests for antibody research, so it is advisable that additional tests be performed for higher sensitivity and for comparison of the results. WAMA Diagnostics manufactures **Imuno-HAI CHAGAS** for indirect hemagglutination and **Imuno-ELISA CHAGAS**.

PERFORMANCE CHARACTERISTICS
Sensitivity: It was evaluated 158 positive specimens. Among 261 negative specimens, it was not found any false positive. Only 5 specimens showed an unspecific fluorescence. This confers to Imuno-Con from WAMA Diagnostics 100% of specificity.

PRECAUTIONS AND WARNINGS
In vitro diagnostic use only.
The kit should be stored in refrigerator at 2-8°C.
Seal the cover slip edge with nail varnish and store in refrigerator at 2-8°C in a dark room for 2 months.
Sodium Azide 0.095% is used and it can react with lead and copper plumbing and highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide buildup. Sodium azide is toxic if ingested.
Human components have been tested for HBsAg, anti-HIV, anti-HCV and are negative. However, this does not assure the absence of these and other infectious diseases. All human components should be handled as potentially hazardous.

Do not interchange kit components with those from sources other than the code number from WAMA DIAGNOSTICA.
Do not use components beyond their expiration date.
Always use in accordance with local regulations.

Periodical maintenance of the microscope is recommended since expired lenses may interfere in the result analysis.
Follow the good laboratory practices (GLP) related to storage, dispensing and material disposal.

DISCONTINUATION
This diagnostic kit replaces this kit since it is not beyond its expiration date. The discontinued kit must be evaluated by Wama's technical support. The discontinued kit will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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