Imuno-CON FTA-Abs Sífilis

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ENGLISH

SUMMARY

Syphilis is a sexually transmitted disease (STD) caused by the spirochete, *Treponema pallidum*.

Other transmission ways are infected blood transfusion, which is rare today due to serological screening, and perinatal (congenital syphilis) that is transmitted

from mother to child in utero.

Primary syphilis is manifested after an incubation period of 10-90 days. During the initial incubation period, patients are asymptomatic. The sore, called a chancre, is a firm, painless skin ulceration occurring in 85% of patients.

Secondary syphilis is characterized after 4 to 10 weeks after primary infection. The most common symptoms are weight loss, headache, anorexia, myalgia, arthralgia, malaise, low fever, lymph nodes and exanthem - occurring in 75 to 100% of patients. On this phase, central nervous system can also be affected. Latent syphilis is after the primary and secondary stages and is characterized by positive serological tests and the absence of clinical findings. This stage can last 1 to 2 years. Tertiary syphilis is manifested in one third of patients without treatment and can present gumma(15%), cardiovascular syphilis (10%) or neurosyphilis (8 - 10%). Serological tests for syphilis are classified as non-treponemal such as VDRL (Venereal Disease Research Laboratory) and RPR (Rapid Plasma Reagin) and treponemal tests such as TPHA (Treponema pallidum Hemagglutination), FTA-Abs (Fluorescent Treponemal Antibody Absorption) and ELISA (Enzyme-Linked Immunosorbent Assay).

FTA-Abs is one of the most sensitive and specific test used on syphilis serology. However, false positive reactions are present in cirrhosis, genital herpes, SLE and, rarely, in pregnant women.

PRINCIPLE OF THE METHOD

Anti-treponemal antibodies on serum bind to an antigen-fixed slide and are identified by fluorescein-labelled, anti-human IgG conjugate.

KIT PRESENTATION

REF 15060-I (60 determinations)

1. 5-well-slide *Treponema pallidum* (12 slides)
2. Fluorescein-labelled, anti-human IgG conjugate (3ml)
3. Phosphate Buffered Saline (PBS) (2 x 50ml)
4. Mounting Medium (4ml)
5. Absorbent solution (Treponema strain Reiter) (3ml)
6. Positive control serum (1ml)
7. Negative control serum (1ml)
8. Coverslips (12 units)
9. Instructions for use

REF 15100-I (100 determinations)

1. 5-well-slide *Treponema pallidum* (20 slides)
2. Fluorescein-labelled, anti-human IgG conjugate (5ml)
3. Phosphate Buffered Saline (PBS) (3 x 50ml)
4. Mounting Medium (4ml)
5. Absorbent solution (Treponema strain Reiter) (5ml)
6. Positive control serum (1ml)
7. Negative control serum (1ml)
8. Coverslips (20 units)
9. Instructions for use

STABILITY AND STORAGE

• **TREPONEMA PALLIDUM SLIDE (1):** Allow the pouch to equilibrate to room temperature for 15 minutes, then remove slide from the pouch. Stable if stored at 2°C to 8°C up to expiration date.

• **FLUORESCIN-LABELED, ANTI-HUMAN IgG CONJUGATE(2):** Ready for use. Stable if stored at 2-8°C up to expiration date. Protect from light. Contains Sodium azide 0.095%

• **PHOSPHATE BUFFERED SALINE (PBS) (3):** Dilute PBS (0.01M; pH7,2 ± 0,2) 1:20 concentrated (Ex: 10ml of PBS + 190ml of distilled or deionized water). Stable if stored at 2-8°C up to expiration date in a clean and sealed container. Discard if turbidity or pH changing is present. Obs.: In case of PBS crystallization before diluting, keep it in water bath at 37°C up to a complete crystal dissolution.

• **MOUNTING MEDIUM (4):** Ready for use. Stable if stored at 2-8°C up to expiration date.

• **ABSORBENT SOLUTION (5):** Intended to remove non-specific anti-treponemal antibody. Ready for use. Shake before using it. Unknown serum should be diluted 1:5 in absorbent solution (Suggestion: 100l of serum + 40ul of absorbent solution). Stable up to expiration date.

• **POSITIVE SERUM CONTROL (6):** Ready for use. Stable up to expiration date. Contains 0.095% Sodium azide.

• **NEGATIVE SERUM CONTROL (7):** Ready for use. Stable up to expiration date. Contains 0.095% Sodium azide.

Obs.: Kit presents results after its first handling and it is stable up to expiration date if stored at 2-8°C.

SPECIMEN COLLECTION AND STORAGE

Only serum specimens should be used for this procedure. Grossly haemolysed, lipaemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store the serum at 2-8°C for up to 72 hours. For long time storage, the serum should be stored at -20°C. Avoid repeated thawing and freezing of samples.

MATERIAL REQUIRED BUT NOT PROVIDED

- Fluorescence microscope and test tubes
- Serological pipette
- Coplin jar or similar
- Tubes 13x75mm and rack
- Distilled or deionized water
- 1 liter container
- Incubation chamber
- Paper towels

PROCEDURE

1. Inactivate the serum for 30 minutes at 56°C.

2. Dilute the unknown serum 1:5 with absorbent solution (See STABILITY AND STORAGE). Allow for 15 minutes to equilibrate to room temperature.

3. Allow the pouche to equilibrate to room temperature for 15 minutes, then remove slide from the pouche. Remove carefully the slide without touching on the substrate. Label it and store into moisture chamber.

4. Add 1 drop of positive (6) and negative controls (7), on wells 1 and 2 respectively and 1 drop of unknown serum on the remaining wells.. Avoid overflowing the wells.

5. Incubate the slide for 30 minutes at room temperature.

6. Remove the slides from the incubation chamber. Hold the slide at tab end and rinse it with 10ml (PBS) (3). Using a pipette, lead PBS by the longitudinal edge of the slide. Be sure not to touch on the reagent wells to avoid damaging substrate. Transfer the slide into Coplin jar or similar and rinse three times with PBS for 5 minutes each time. Gently shake the jar during rinsing.

7. Remove the slide from Coplin jar, one by one. Blot the edge of slide on a paper towel to remove PBS. Dry around reactive wells. To prevent slide from drying, go IMMEDIATELY to step 8.

8. Place the slide back to incubation chamber. Apply 1 drop of Fluorescein-labeled, anti-human IgG conjugate on each well. Be sure to fill the well completely.

9. Incubate the slide for 30 minutes at room temperature. Protect it from light.

10. Repeat step 6.

11. Remove the slide from Coplin jar, one by one. Blot the edge of the slide on a paper towel to remove PBS excess. Dry around the reactive wells. To prevent wells from drying, go IMMEDIATELY to step 12.

12. Apply 3 to 4 drops of Mounting Medium (4) on reactive wells. Cover the slips with coverslips avoiding air bubbles. Dry the Mounting Medium excess with a paper towel. Clean the slide reverse.

13. The slide must be read under a fluorescence microscope as soon as possible. Otherwise, store the slide at 2-8°C, protect it from light and read it in the following day. The Mounting Medium should not dry, if necessary, apply more.

INTERPRETATION OF RESULTS

NEGATIVE: ABSENCE of yellow-greenish fluorescence on spirochete (*T. pallidum*). Spirochete is not easily seen.

Positive: PRESENCE of a typical yellow-greenish fluorescence on spirochete (*T. pallidum*).

OBSERVATION: Always check positive and negative controls for reaction control purposes.

INTERPRETATION

Negative serum are considered as **NEGATIVE**, and they are represented by non infected patients with titers lower than 1/5. Positive serum are considered **POSITIVE** and they are represented by infected patients by *Treponema pallidum* or previous infection with titers = or > to 1/5. The majority of the cases, even after the disappearance of the disease, FTA-Abs is positive for some years.

Obs.: It is not recommended to submit positive samples to titration since through this method since high titers can persist for several years. This does not mean the patient is in any infection phase, that is, other methods such as VDRL or RPR can indicate better treatments associated to the patient's health status. If the titration is needed, use indirect hemagglutination method (TPHA).

METHOD LIMITATION

The use of this method avoid possible unpecific answers, therefore, the results should be used in conjunction with available information from clinical evaluation and other diagnostic procedures.

The source of light, filters and optic from different microscopes of fluorescence can affect the test sensitivity.

PERFORMANCE CHARACTERISTICS

Clinical and diagnostic sensitivity: 100% of sensitivity. 27 in-house positive specimens were compared and confirmed. All results were satisfactory. No false negative was found.

Specificity: 100% of specificity. 273 in-house negative specimens were compared and confirmed. All results were satisfactory. No false positive was found.

PRECAUTIONS AND WARNINGS

1. For in vitro diagnostic use only.
2. The kit should be stored at 2°C to 8°C.
3. For a better slide storage, seal the coverslip edge with nail varnish and store at 2-8°C in a dark room up to 2 months.
4. Sodium Azide 0.095% is used and it can react with lead and cooper plumbing forming highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide buildup. Sodium azide is toxic if ingested.
5. All human components have been tested for HBsAg, anti-HIV, anti-HCV and found to be negative. However, this does not assure the absence of these and other infectious diseases. All human components should be handled as potentially hazardous.
6. Do not interchange kit components with those from sources other than the same code number from WAMA Diagnostica.
7. Do not use components beyond their expiration date.
8. Disposal in accordance with local regulations.
9. The periodical maintenance of the microscope is recommenden since on expired lamp can interfere in the result analysis.
10. Follow the good laboratory practices (GLP) related to storage, dispensing and material disposal.

WARRANT

WAMA Diagnóstica replaces this kit since it is not beyond expiratin date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replace if technical support finds

evidence that running, houlding and storage were not properly fallow ed.

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