# Imuno-CON ANA-HEp-2

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## ENGLISH

#### SUMMARY

Antinuclear antibodies (ANA) are present in the serum of patients with connective tissue

disorders.

Antinuclear antibodies (ANA) are related to certain disturbs and they are regarded as markers for some diseases. For example, antibodies against Smith antigen (anti-Sm) and against DNA (anti-DNA) are highly associated to Systemic Lupus Erythematosus (SLE) whereas the anti-centromere antibodies are associated to a limited form of scleroderma, variant CREST. ANAs are regarded as auto-antibodies because in some rheumatic diseases they are also directed to components of extranuclear cells such as cytoplasmatic elements and cyto-skeletal.

In 1982, American Rheumatism Association included a positive test for ANA as a criterion for Systemic Lupus Erythematosus (SLE) diagnosis. Therefore, a negative result discard SLE, however, it is not specific, since patients with connective tissue disorders such as arthritis rheumatoid, sclerodermia, Sjogren Syndrome and dermatomiositis are frequently ANA positive. Besides, positive ANA with low titers can also be found in infections, neoplasia and on healthy elderly people.

Indirect immunofluorescence method is considered to be the standard for ANA, in which tissue cut (mouse's liver or kidney) or human epithelial cells (HEp-2) is used as substrate. HEp-2 cells, unlike tissue cut, express almost all antigens, including centromere, Ro (SSA), La (SSB), Scl-70, ANCP etc.

The reagents of this kit are aimed at detecting important antibodies by using HEp-2 cells and IgG conjugate as substrate which eliminate physiological false positive results normally present due to low titers of specific IgM.

**PRINCIPLE OF THE METHOD**
Antinuclear antibodies present in the serum bind to the antigenic substrate (HEp-2 cells) and they are detected by an antigramma globulin labelled with FITC.

#### KIT PRESENTATION

**REF 121004** (100 determinations)

- 10-well HEp-2 cells Slide (10 slides).
- FITC IgG conjugate (1x5ml).
- Phosphate Buffered Saline (PBS) 20x (1x50ml).
- Mounting Medium (1x4ml)
- Homogeneous positive control serum (1x0,5ml)
- Centromere positive control serum (1x0,5ml)
- Negative control serum (1x1ml)
- Evans blue (1x2ml)
- Cover slips (10 unidades)
- Package insert

**REF 122004** (200 determinations)

- 10-well HEp-2 cells Slide (20 slides).
- FITC IgG conjugate (2x5ml).
- Phosphate Buffered Saline (PBS) 20x (2x50ml).
- Mounting Medium (2x4ml)
- Homogeneous positive control serum (1x1ml)
- Centromere positive control serum (1x1ml)
- Negative control serum (1x1ml)
- Evans blue (1x2ml)
- Cover slips (20 unidades)
- Package insert

#### REAGENT STABILITY AND STORAGE

- 10-WELL HEp-2 CELLS SLIDE (1):** Let the pouch reach to room temperature for 15 minutes, then remove the slide. Stable if stored at 2-8°C up to expiration date.
- FITC IGG CONJUGATE (2):** Ready for use. Stable if stored at 2-8°C up to expiration date. Protect from light. It contains sodium azide 0.095%
- PHOSPHATE BUFFERED SALINE (PBS)20x (3):** Dilute PBS on distilled or deionized water. 1 PBS + 19 distilled or deionized water. Stable for 1 month if stored at 2-8°C in a clean and sealed container. Discard if turbidity is present. Diluted PBS is used as washing buffer and diluent for patient specimens.

Stable in refrigerator for 12 months. Ready for use. Stable if stored at 2-8°C up to expiration date. **DO NOT STORE** below 2°C to avoid crystallization. It contains sodium azide 0.095%.

**•HOMOGENEOUS POSITIVE CONTROL SERUM (5):**Ready for use. Store at 2-8°C. Stable up to expiration date. It contains sodium azide 0.095%

**•CENTROMERE POSITIVE CONTROL SERUM (6):** Ready for use. Store at 2-8°C. Stable up to expiration date. It contains sodium azide 0,095%.

**• NEGATIVE CONTROL SERUM (7):** Ready for use. Store at 2-8°C. Stable up to expiration date. It contains sodium azide 0.095%

**• EVANS BLUE (8):** Ready for use. Store at 2-8°C. Stable up to expiration date.

**Note.** The kit retains the same performance after first use and is stable until the expiration dara described on the label, provided it is kept at the indicated temperature (2-8 °C).

#### SPECIMEN COLLECTION AND STORAGE

Only serum specimens should be used for this procedure. Grossly haemolysed, lipaemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store serum at 2-8°C up to 72 hours. For longer storage, the serum should be stored in a freezer.

Avoid repeated thawing and freezing of specimens.

#### MATERIAL REQUIRED BUT NOT PROVIDED

- Fluorescence microscope
- Micropipette
- Test tube
- Serological pipettes
- Coplin jar or similar
- Distilled or deionized water
- 1 liter container
- Incubation chamber
- Papertowel

#### PROCEDURE

**Qualitative Test (screening):** It is intended to screening and removal of negative serum.

- Dilute patient serum 1:40 with diluted PBS (3) (10µl of serum + 390µl of diluted PBS).
- Allow the pouch reach to room temperature for 15 minutes, then remove slide. Carefully remove the slide without touching on the substrate. Label it and store in the moisture chamber.
- Add 1 drop of homogeneous positive control (5), centromere positive (6) and negative control serum (7), undiluted, in the wells 1, 2 and 3 respectively. Avoid overfilling the wells.
- Add 1 drop of unknown serum diluted at 1:40 in the remaining wells. Avoid overfilling the wells.
- Incubate the slide for 30 minutes at room temperature.
- Remove the slide from the incubation chamber. Hold the slide at tab end and rinse it using 10ml PBS (3). With a pipette, label PBS by the longitudinal edge of the slide. Be sure not to touch reagent wells to avoid damaging the substrate. Transfer the slide into the Coplin jar or similar and wash it twice with diluted PBS (3), for 5 minutes, gently shake the jar.
- Remove the slide from the Coplin jar, one by one. Blot the edge of the slide on a paper towel to remove PBS excess. Dry around reactive wells. To prevent slide wells from drying, go

IMMEDIATELY to step 8.  
**SUBSTRATE DRYING MAY AFFECT THE RESULTS.**  
 8. Place the slide back to incubation chamber. Apply 1 drop of FITC IgG conjugate (2) in each well. Be sure to fill the well completely.  
 9. Incubate the slide for 30 minutes at room temperature. Protect it from light.  
 10. Remove the slide from incubation chamber then transfer it into a Coplin jar or similar with PBS to remove conjugate excess. Transfer the slide to another Coplin jar or similar with diluted PBS and wash twice for 5 minutes each. Gently shake the jar. Apply 2 to 3 drops of Evans blue (8) into Coplinjar in this last washing.  
 11. Remove the slide from Coplin jar, one by one. Blot the edge of the slide on a paper towel to remove PBS excess. Dry around reactive wells. To prevent wells from drying, go IMMEDIATELY to step 12.  
 12. Apply 3 to 4 drops of Mounting Medium (4) in the reactive wells. Cover the slide with cover slip. Avoid air bubbles. Dry the excess of Mounting Medium with a paper towel.  
 13. Slide must be read under a fluorescence microscope as soon as possible. However, if it is not possible, store the slide at 2-8°C and protect it from light (dark room) for up to a week.

b. **Semi-quantitative Test.** (Titration): It is intended to determine the titers of antibodies for positive serum through the screening test.  
 1. Dilute positive serum with diluted PBS (3) following ratio 2. Ex: 1:40, 1:80, 1:160, 1:320, ... 1:2560 and so on.  
 2. Allow the pouch equilibrate room temperature for 15 minutes, then remove the slide. Carefully remove the slide without touching on the substrate. Label it and store in the incubation chamber.  
 3. Add 1 drop of homogeneous positive control (5), centromere positive (6) and negative (7), undiluted, on1, 2 and 3 wells respectively. Avoid overfilling the wells.  
 4. Add 1 drop of patient's diluted serum in the other wells, using one well for each dilution. Avoid overfilling the wells.  
 5. Follow the steps 5 through 13 of the screening tests.

**READING OF RESULTS**

**Negative:** ABSENCE of a typical apple-green fluorescence color. The nuclear image of the specimen should be always similar to that of negative control.

**Positive:** PRESENCE of a typical apple-green fluorescence color. The nuclear and/or cytoplasmatic image of the specimen is more fluorescent than that of negative control.

**INTERPRETATION.**

**1. Homogeneous pattern:** A solid color of all nucleus in interphase, with occasional negative nucleus. Chromosomal region of the cells in metaphase is positive. Nuclear antigens in decrescent order of frequency are: DNP, histone, DNA double/simple strip (dsDNA/ssDNA). **Interpretation:** High titers are suggestive of idiopathic SLE and a test for anti-DNA more specific and quantitative should be performed. Low titers are present, occasionally, in rheumatoid arthritis, idiopathic LES and induced by drugs.

**Peripheral pattern:** A solid color around the nuclear membrane with weak color in the core of the nucleus, in interphase. Chromosomal region of mitotic cells is positive, like homogeneous pattern. Antigens are DNA double spiral and rarely DNA simple strip, DNP or histone.

**Interpretation:** High titers are suggestive of idiopathic LES. Specific and quantitative determinations of anti-DNA should be performed. Low titers suggest other conjunctive tissue diseases.

**3. Dotted pattern:** Granulou or thin color, uniform or variable size of the nucleus in interphase is visible, normally colorless nucleolus. The chromosomes are negative, but the outside region of the chromosomal area presents variable color depending of antigens. Nuclear antigens are Sm, RNP, Scl-70, SS-B/La, SS-A/Ro and many other antigen-antibody systems.

**Interpretation:** Sm presence suggests SLE. RNP may be present in combinated diseases of the conjunctive tissue, myositis, scleroderma or SLE, Scl-70 in the scleroderma, SS-B/La and SS-A/Ro in the Sjogren syndrome, although SS-A/Ro can present different patterns of color.

**4. Nucleolar pattern:** Thick granulous inside the nucleus is visible, normally less than 6 per cell with or without thin granulous. The chromosomal region of mitotic cells is negative. Nuclear antigens can present different sites with cellular cycle. Antigens are Pm-Scl, nucleolin, Ku and others unknown.

**Interpretation:** It occurs on Reynaud phenomenon and other ways of scleroderma and its variant, but rarely on SLE. High titers can be linked, besides scleroderma, to Sjogren syndrome.

**5. Centromere pattern:** Only pattern regarded as defined. A slight granular color, multiple of 46, very discreet and observed. In the Hep-2 occurs a typical color of mitotic cells in the cromosomal area and in the interphase nucleus.

**Interpretation:** Anti-centromere antibody is highly suggestive of a scleroderma variant known as CREST (calcinosis, Reynaud fenomeno, dissimilitiy of oesophagus, sclerodactilia and telangiectasis)

**Mitochondrial pattern:** Irregular thick granules, in filaments format through all cytoplasm. Antigen is M2.

**Interpretation:** This pattern is typical in primary bilial cirrose and rare on conjunctive tissue disease.

**METHOD LIMITATION**

In some cases, positive serum for nuclear fluorescence are weak reagent or negative in the initial dilution (prozone effect). In such doubtful cases, the serum should be screened to higher dilutions, and, if positive, the titers of the antibodies should be determined. A high titer for ANA suggests a conjunctive tissue disease. But not all positive ANA should be regarded as diagnostic. The test results should be evaluated together with other serological results as well as patient's clinical assay. The same specimen can present two or more types of fluorescence patterns. This phenomenon can be observed when the specimen is stored up to its end point.

The incidence of positive ANA increases with the patient's age. This incidence can also occur in a small proportion of patients with infectious and/or neoplastic diseases.

Combinations of auto-antibodies can induce to homogeneous and/or dotted patterns and therefore, it is recommended that specific tests for anti-dsDNA and ENA are performed with all homogeneous pattern specimens.  
 Fluorescence microscopes with different light source, filters and optic and also from different brands can affect sensitivity of the tests.

**PRECAUTIONS AND WARNINGS**

1. For in vitro Diagnostic use only.  
 2. The kit should be stored at 2-8°C.  
 3. The conservation of the slide is done by sealing the coverslip edge with nail varnish and storing at 2-8°C in a dark room, for 2 months.  
 4. Sodium Azide 0.095% is used and it can react with lead and cooper plumbing forming highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide buildup. Sodium azide is toxic if ingested.

5. All human components have been tested for HBSAg, anti-HIV, anti-HCV and found to be negative. However, this does not assure the absence of these and other infectious diseases. All human components should be handled as potentially hazardous.

6. Do not interchange kit components with those from sources other than the same code number from WAMA/DIAGNOSTICA.  
 7. Do not use components beyond their expiration date.  
 8. Dispose in accordance with local regulations.  
 9. The periodical maintenance of the microscope is recommended since an expired lamp can interfere in the result analysis.  
 10. Follow good laboratory practices (GLP) related to storage, dispensing and material disposal.

**WARRANT**  
 WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replace if technical support finds evidence that running, handling and storage were not properly followed.

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- Ben-Chetrit, E.: The molecular basis of the SS-A/Ro antigens and the clinical significance of their autoantibodies. *Br. J. Rheum.* 32: 396-402, 1993.
- Fritzer, M. J.; Kinsella, T. D. and Garbutt, E.: The CREST syndrome: A distinct serologic entity with anti-centromere antibodies. *Am. J. Med.* 69: 520, 1980.
- Humbel, R. L.: Detection of antinuclear antibodies by immunofluorescence. In: Van Venrooij, W. J. ; Maini, R. N. (editor): **Manual of Biological Markers of Disease.** Kluwer Academic Publisher, A21, A21/6, 1993.
- Molden, D. R.; Nakamura, R. M. and Tan, E. M.: Standardization of immunofluorescence test for autoantibody to nuclear antigens (ANA). Use of reference sera of defined antibody specificity. *Am. J. Clin. Pathol.*, 87: 57, 1984.
- Naparstek, Y.: The role of autoantibodies in autoimmune disease. *Ann. Rev. Immunol.*, 11: 79, 1993.
- Tan, E. M.; Cohen, A. S.; Fries, J. F.; et al.: The 1982 revised criteria for the classification of systemic lupus erythematosus. *Arthritis Rheum.*, 25(11), 1271-77, 1982.
- Tan, E. M.: Autoantibodies to nuclear antigens (ANA). Their immunobiology and medicine. *Adv. Immunol.*, 33: 167-239, 1982.
- Tan, E. M.: Antinuclear antibodies: diagnostic markers for autoimmune for diseases and probes for cell biology. *Adv. Immunol.*, 40: 93-151, 1989.
- Van Venrooij, W. J.; et al.: The consensus workshops for the detection os autoantibodies to intracellular antigens in rheumatic diseases. *J. Imm. Meth.*, 140: 181-189, 1991.
- Van Muihlen, C. A. and Tan, E. M.: Autoantibodies in the diagnosis of systemic rheumatic diseases. *Sem. Arthr. Rheum.*, 24(5): 323-358, 1995.

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