Imuno-Elisa CHAGAS

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ENGLISH

SUMMARY

Chagas disease or American Trypanosomiasis is an endemic and chronic disease. Its pathogenic agent is a flagellate protozoan called *Trypanosoma cruzi* which is transmitted to humans by a bug, triatominae.

3stages can be present: acute, latent or indetermined and chronic.

The acute phase can present manifestations that include fever, myocardiopathy,

lymphadenopathy, hepatosplenomegaly and parasitaemia.

Pseudocyst is formed after intracellular multiplication of parasite on muscles and on cells of endothelial reticulum system.

Latent and intermediary stages are asymptomatic. The disease can evolve for the chronic phase with manifestations of myocardiopathy, degeneration of ganglion cells of peripheral and central nervous system and hypertrophy and dilatation of some organs such as oesophagus and colon.

Chagas disease remains a major health problem in many Latin American countries because the high numbers of deaths.

As the minority of patients presenting positive sorology for *T. cruzi* develop clinical evidence for chronic disease, information provided by the laboratory has a decisive and fundamental role for etiologic diagnosis.

Chagas disease methods are: Complement fixation, agglutination, precipitation, immunofluorescence, hemagglutination and ELISA from which the most used are indirect hemagglutination (HAI), indirect immunofluorescence (IFI) and ELISA.

ELISA methodology is increasingly used due to its easy-to-use format, excellent sensitivity and specificity and the possibility of automation.

PRINCIPLE OF THE METHOD

The microtiter wells of the plate are covered with highly purified recombinant antigens of *Trypanosoma cruzi*.

The specific antibodies Anti-*Trypanosoma cruzi* present in serum bind to these recombinant antigens. The unbound material is removed by washing and an antigammaglobulin anti-IgG human marked with peroxidase is added to the reaction. This conjugate binds to the specific antibodies anti-*Trypanosoma cruzi* IgG human binded on the antigen on the plate.

The unbound material is removed again by washing and a substrate (TMB) is added which will develop a color on the wells in the presence of enzyme (peroxidase), indicating the presence of human antibody anti-*Trypanosoma cruzi*. The enzymatic reaction is stopped by the addition of a stop solution and the absorbance is measured at 450 nm. The concentration of specific IgG antibody is directly proportional to the intensity of the color of the reaction.

KIT PRESENTATION

REF 44096-E (96 determinations)

- Microplate with antigen of *Trypanosoma cruzi* (12x8 wells).
- Diluent solution "orange" (1 x 100ml)
- Wash buffer - 20 x concentrated (1 x 50 ml)
- Conjugate anti-human IgG Peroxidase (11ml)
- Cromogen (TMB) (11ml)
- Stop solution (11ml)
- Negative serum control (2ml)
- Positive serum control low (2,5ml)
- Positive serum control high (2ml)

10. Instructions for use

REF 44192-E (192 determinations)

- Microplate with antigen of *Trypanosoma cruzi* (12x8 wells).
- Diluent solution "orange" (2 x 100ml)
- Wash buffer - 20 x concentrated (2 x 50ml)
- Conjugate anti-human IgG Peroxidase (2 x 11ml)
- Cromogen (TMB) (2 x 11ml)
- Stop solution (2 x 11ml)
- Negative serum control (2 x 2ml)
- Positive serum control low (2 x 2,5ml)
- Positive serum control high (2 x 2ml)

10. Instructions for use

MATERIAL REQUIRED BUT NOT PROVIDED

- Micropipettes
- Distilled or deionised water
- Incubator at 37°C.
- ELISA plate reader with 450nm filter
- Paper towel.
- Immunological container.

STABILITY AND STORAGE

MICROPLATE (1): Stable if stored at 2-8°C up to the expiration date. Remove from the envelope the number of strips needed and let them reach to room temperature (20-25°C). Keep the remaining strips in the aluminium pouch with desecant at 2-8°C.

DILUENT SOLUTION (2): Stable if stored at 2-8°C up to the expiration date. This solution is should be tested within 8 hours. Do not freeze. It contains *Proclin 300-0.05% as preservative. Allow the components to reach room temperature (20-25°) prior to use.

WASH BUFFER (3): Dilute wash buffer using 1 part of the buffer with 19 parts of distilled or deionised water. During each washing cycle, each well should be filled with 300µl. Wash buffer is diluted and stable for 4 weeks at 2-8°C. It contains *Proclin 300-0.05% as preservative. Allow the components to reach room temperature (20-25°) prior to use.

CONJUGATE ANTI-HUMAN IgG PEROXIDASE (4): Ready for use. Stable if stored at 2-8°C up to the expiration date. Do not freeze. It contains *Proclin 300-0.05% as preservative. Allow the components to reach room temperature (20-25°) prior to use.

CROMOGEN (TMB) (5): Ready for use. Stable if stored at 2-8°C up to expiration date. Do not freeze. It contains *Proclin 300-0.05% as preservative. Allow the components to reach room temperature (20-25°) prior to use. TMB is not mutagenic and carcinogenic.

STOP SOLUTION (6): Ready for use. This solution is hydrochloric acid (Hcl), therefore, it is corrosive to check carefully. Wash with abundant water. Allow the components to reach room temperature (20-25°) prior to use. Stable up to the expiration date.

NEGATIVE SERUM CONTROL (7): Ready for use. Do not dilute. Store at 2-8°C up to expiration date. It contains *Proclin 300 - 0.05% as preservative.

POSITIVE SERUM CONTROL LOW (8): Ready for use. Do not dilute. Store at 2-8°C up to the expiration date. It contains *Proclin 300 - 0,05% as preservative. It must be used in duplicate to check the best performance and cut-off value.

POSITIVE SERUM CONTROL HIGH (9): Ready for use. Store at 2-8°C up to the expiration date. It contains *Proclin 300 - 0.05% as preservative.

* Proclin 300 is ROHM&HAAS Limited registered trademark.

The kit presents good performance after being used for the first time. It is stable up to the expiration date if store at 2 - 8°C.

SPECIMEN COLLECTION AND HANDLING

Serum specimens free of haemolysis, lipemia and contamination should be used. Plasma specimens with EDTA, heparin or oxalate can interfere on the results and should be avoided.

Serum can be stored at 2-8°C for 48 hours. For longer storage, keep at -20°C up to 1 year.

Frozen samples must be homogenised prior to test. Avoid air bubbles. Do not use Sodium azide as preservative since the enzyme peroxidase can be inhibited.

After diluting the serum with diluent solution, use it within 8 hours. Avoid repeated freezing and thawing cycles.

DILUTION OF SPECIMEN

Dilute the serum to 1:25.

Tube

20µl of serum + 480µl of diluent solution. Best within 8 hours. Dispense 125µl of diluted specimens into the plate wells.

Plate

5µl of serum + 120µl of diluent solution and proceed the test as per described.

CONTROL

Controls should be included to guarantee precise interpretation of the results. The controls are ready for use. Do not dilute. Dispense them into the plate wells as per distribution plan.

PROCEDURE

1. According to the established distribution plan, use 1 well for negative serum (7), 1 well for positive control high(9) and 2 wells for positive control low(8), which will determine the cut-off value.

2. Dispense 125µl of controls and diluted specimens at 1:25 in case of dilution in tube. With the fingers, beat on the plate edges or by mechanical vibration for 15 seconds for an homogeneous mixture.

3. Cover the plate with an adhesive sheet, place the plate on a humid paper towel and incubate at 37°C for 60 minutes.

4. Discard the material in a container with disinfectant solution.

5. Wash three times the plate with diluted wash buffer (3) (see Stability and Storage steps). Fill and empty the well with 300µl of wash buffer. A manual or automatic microplate washer is recommended.

6. After the last washing, invert the plate and beat it against an absorbent paper to remove residual liquids.

7. Dispense 100µl of conjugate (4) into each well. With fingers, beat on the plate edges or by mechanical vibration (plate agitator) for 15 seconds for an homogeneous mixture.

8. Cover the plate with an adhesive sheet, place the plate on a humid absorbent paper and incubate at 37°C for 30 minutes.

9. Repeat the steps 4, 5 and 6.

10. Dispense 100µl of cromogen (5) onto each well. With fingers, beat on the plate edges or by mechanical vibration (plate agitator) for 15 seconds for an homogeneous mixture.

11. Incubate the plate at 37°C on a dry surface for 15 minutes. Protect from light.

12. Block the reaction by dispensing 100µl of stop solution (6) into each well. With the fingers, beat on the well edges or by mechanical agitation (plate agitator) for 30 seconds for an homogeneous mixture assuring that the blue color change to yellow.

13. After blocking of reaction, read IMMEDIATELY using a microplate reader with 450nm filter.

WASHING:

• Washing step is essential for the test performance. Insufficient washings will result in imprecised readings or absorbances falsely high.

• It is recommended to duplicate the tests.

• Do not take more than 30 minutes to dispense the specimens and controls.

SUMMARY OF PROCEDURE

1. Dilute the serum 1:25 (20µl of serum + 480µl of diluent solution (2));

2. Dispense 125µl of diluted serum and controls into each well;

3. Mix for 15 minutes and incubate at 37°C for 60 minutes;

4. Discard the plate content and wash 3 times with wash buffer (3);

5. Dispense 100µl of conjugate (4) into each well;

6. Mix for 15 seconds and incubate at 37°C for 30 minutes;

7. Discard the plate content and wash for 3 times with wash buffer (3);

8. Dispense 100µl of cromogen (5) into each well.

9. Mix for 15 seconds and incubate at 37°C for 15 minutes in a dark room.

10. Dispense 100µl of stop solution (6) into each well.

11. Mix for 30 seconds;

12. Read OD using a microplate reader with 450nm filter.

RESULTS

Determine OD for each test and serum control.

Cut off = OD Mean of positive serum control low

Validation of Test: OD Mean of positive serum control low should be 3 times higher than the OD of Negative Serum control.

INTERPRETATION

POSITIVE : OD of specimen is higher than borderline (dubious zone).

DUBIOUS RESULT (Borderline): OD of specimen is lower than O.D. of positive control serum low, but is higher than Cut-off value. It is called the dubious zone.

NEGATIVE : OD of specimen is lower than Cut-off value.

TEST PERFORMANCE

It was made tests to determine the sensitivity and specificity of Imuno-ELISA CHAGAS from WAMA Diagnostica. 120 known positive specimens were used and any false negative result was found. It awarded to Imuno-ELISA CHAGAS from WAMA Diagnostica a sensitivity of 100% on the specificity of the test. It was used 642 known negative specimens and any positive result was found what gave a specificity of 100% to Imuno-ELISA from WAMA Diagnostica.

LIMITATION OF USE

For a correct interpretation, patient's clinical data should be taken in account.

It is recommended the use of other methodology for the Chagas interpretation results. WAMA manufactures ImunoCon CHAGAS (Code. 1460-I, 14100-I, 14200-I) indirect immunofluorescence method and Imuno-HAI CHAGAS (Code. 34096-H, 34192-H, 34380-H) indirect hemagglutination method.

PRECAUTIONS AND WARNINGS

1. Reagents should be stored at 2-8°C

2. The expiration date is the last day of the month printed on the labels and on the kit box.

3. Avoid high temperatures and sunlight exposition.

4. Do not freeze the reagents. Irreversible damage.

5. Reagents contain Proclin 300 (0.05%) as preservative. Toxic if ingested. In case of skin contact, wash with plenty quantities of water.

6. The source material of Imuno-ELISA CHAGAS from WAMA used in the preparation of the controls

has been tested and confirmed negative for HCV, HIV and HBsAg. Because no test can offer complete assurance that products derived from human source will not transmit infectious agents it is recommended that these controls be treated as potential biohazards materials.

7. Do not interchange reagents or microplates from different lots.

8. Allow the reagents to reach room temperature (20-25°C) prior to use.

9. Reagents should be homogenized by inversion and sealed immediately after use.

10. Avoid the drying of wells during the test.

11. Avoid contamination by repeatedly pipetting the reagents.

12. Samples and controls should be tested and used at the same time to maintain test conditions.

13. Do not use components after the expiration date.

14. Dispose in accordance with the local regulations.

15. Follow good manufacturing practices (GMP) at storing, handling and material disposable.

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