Imuno-HAI CHAGAS

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ENGLISH

SUMMARY

Chagas disease or American Trypanosomiasis is an endemic and chronic disease. Its pathogenic agent is a flagellate protozoan called *Trypanosoma cruzi* which is transmitted to humans by a bug, triatominae. 3 stages can be present: acute, latent or indetermined and chronic.

The acute phase can present manifestations that include fever, myocardiopathy, lymphadenopathy, hepatosplenomegaly and parasitaemia.

Pseudo-cystus is formed after intracellular multiplication of parasite on muscles and on cells of the endothelial reticulum system.

Latent and intermediary stages are asymptomatic. The disease can evolve for the chronic phase with manifestations of myocardiopathy, degeneration of ganglion cells of peripheral and central nervous system, hypertrophy and dilatation of some organs such as esophagus and colon.

Chagas disease remains a major health problem in many Latin American countries because the high numbers of deaths.

As the minority of patients presenting positive sorology for *T. cruzi* patients develop clinical evidence for chronic disease, the information provided by the laboratory has a decisive and fundamental role for etiologic diagnostic.

Chagas disease methods are: Complement fixation, agglutination, precipitation, immunofluorescence, hemagglutination and ELISA from which the most used are indirect hemagglutination (HAI), indirect immunofluorescence (IFI) and ELISA.

PRINCIPLE OF THE METHOD

Trypanosoma cruzi sensitised formalised tanned fowl erythrocytes. These erythrocytes reveal agglutination when react with antibodies against these antigens present in patients serum.

KIT PRESENTATION

REF 34096-H (96 qualitative tests)

1. *Trypanosoma cruzi* sensitised erythrocytes (1 x 2.4 ml)
2. Diluent solution (1 x 40ml)
3. 2-Mercaptoethanol (0,5ml)
4. Positive serum control (1ml)
5. Negative serum control (1ml)
6. Disposable "V"-well microtitration plate (1x96 wells)
7. Instructions for use

REF 34192-H (192 qualitative tests)

1. *Trypanosoma cruzi* sensitised erythrocytes (1 x 4.8 ml)
2. Diluent solution (1 x 60ml)
3. 2-Mercaptoethanol (0,5ml)
4. Positive control serum (1ml)
5. Negative control serum (1ml)
6. Disposable "V"-well microtitration plate (2 x 96 wells)
7. Instructions for use.

REF 34380-H (380 qualitative tests)

1. *Trypanosoma cruzi* sensitised erythrocytes (2 x 4.8 ml)
2. Diluent solution (2 x 60ml)
3. 2-Mercaptoethanol (1ml)
4. Positive control serum (1ml)
5. Negative control serum (1ml)
6. Disposable "V"-well microtitration plate (4 x 96 wells)
7. Instructions for use.

REAGENT STABILITY AND STORAGE

• **TRYPANOSOMA CRUZI SENSITISED ERYTHROCYTES (1):** Ready for use. Stable if stored at 2-8°C up to expiration date. **Do not freezer.** Homogenize before using, always in circular movements and for vial inversion but never agitate, since this movement can induce to haemolysis and/or agglutination. Allow the components reach to room temperature prior to use.

• **DILUTION SOLUTION (2):** Keep at 2-30°C. When ready to test, add 70 µl of 2-mercaptoethanol for 10ml of diluent solution. This mixing is intended to dilute the serum specimen . The use of 2-mercaptoethanol is intended to remove the antibodies which cause interference.

• **POSITIVE CONTROL SERUM (4):** Ready for use. Stable if stored at 2-8°C up to expiration date.

• **NEGATIVE CONTROL SERUM (5):** Ready for use. Stable if stored at 2-8°C up to expiration date.

Obs.: Kit presents results after its first handling and it is stable up to expiration date if stored at 2-8°C.

SPECIMEN COLLECTION AND STORAGE

Fresh serum from fasting patients or stored at -20° should be used. Do not repeatedly freeze-thaw the serum. Old serum can freeze in contact with 2-mercaptoethanol, causing false positive results. Inactivated serum is recommended.

MATERIAL REQUIRED BUT NOT PROVIDED

- Test tubes
- Serological pipettes
- Glasswire
- Disposable material container

PROCEDURE

A. Qualitative test

1. Lay the plate onto a humid piece of cloth to avoid electrostatics.
 2. Use 1 well of the plate for each specimen, including positive and negative controls.
 3. Dilute the serum at 1:32 together with diluent solution (2) and 2-mercaptoethanol (3). Suggestion: Dilute 10 µl of specimen + 0.31ml of diluent solution with 2-mercaptoethanol in a test tube.
- Obs.:** Positive and negative control serum are ready for use (pré-diluted at 1:32). They should not be diluted.
4. Pipette 50 µl of positive and negative controls and serum diluted at 1:32 of each sample on the respective wells. Suggestion: At 1 positive control, A2

negative control and A3, A4, A5...sera to be tested.
5. Add 25µl of erythrocytes homogeneous suspension (1) in each well.
6. Shake the plate by mechanical vibration (plate agitator) or beating with the fingers on the edges of the plate for 3 to 4 minutes.
7. Let rest for 1 to 2 hours at room temperature in a place free of vibration (IMPORTANT).
8. Read the results.

READING

NEGATIVE: Erythrocytes settle in the bottom of the well forming a button.
Positive: Erythrocytes settle in the bottom of the well, like a rug, often presenting irregular edges.

B. Semi-quantitative test

The serum which shows positive reaction in the qualitative test should be titered.

1. Dilute the serum at 1:32 according to step 3 of qualitative test.
2. Pipette 50µl of diluent solution with 2-mercaptoethanol (3) from the second well of the plate up to dilution that is intended to evaluate. Ex: If the dilution of the serum is up to 1:512, pipette 50µl of diluent in the wells A2, A3, A4 and A5. Do not pipette in the well 1.
3. Transfer 50µl of serum diluted at 1:32 in the wells A1 and A2. Homogenize the serum with the diluent in the well A2 (dilution 1:64) and transfer 50µl in the well A3 (dilution 1:128) and so on, discard 50µl in the end.
4. Pipette 25µl of erythrocytes homogeneous suspension (1) in all wells.
5. Shake the plate by mechanical vibration (plate agitator) or beating with the fingers on the edges of the plate for 3 to 4 minutes.
6. Let rest for 1 to 2 hours at room temperature in a place free of vibration. (IMPORTANT).
7. Read the results.

INTERPRETATION

Titers < 1:32: NEGATIVE

Titers = or > 1:32: POSITIVE

TITER OF THE SPECIMEN: The highest dilution showing positive reaction. The end point is when the erythrocytes cover 50% of the bottom of the well.

NOTE: It is advisable to use more than one kind of test in order to have a better sensitivity and results comparison when researching anti-*T.cruzi* antibodies. WAMA manufactures the kit **Imuno-CON CHAGAS**, immunofluorescence method and **Imuno-ELISA CHAGAS**, ELISA method.

PERFORMANCE CHARACTERISTICS

Clinical and diagnostic sensitivity: 100% of sensitivity, 83 in-house positive specimens were compared with Hemacruzi kit - BioMerieux. All results were satisfactory. No false negative was found.

Specificity: 98.7 of specificity, 235 in-house negative specimens were compared with Hemacruzi kit - BioMerieux. All results were satisfactory. It was found only 4 false positive results.

Precision

Intra-assay: The precision was verified by using 10 specimens in duplicate and performance in 3 different days. The negative and positive values were identified at 100% of the times.

Inter-assay: The precision was verified by using the same 10 specimens and performance 3 assays with 3 different lots under the same conditions (date and operator). The negative and positive values were identified at 100% of the times.

PRECAUTIONS AND WARNINGS

1. Do not use reagents from different lots.
2. The plates should be free of vibration since it can affect the interpretation of results
3. The vials of the erythrocytes should be let stand to avoid irreversible auto-agglutination.
4. USE EXCLUSIVELY THE PLATE OF THE KIT. It should not be reused.
5. Cross reaction with visceral and legumentar Leishmaniasis can be present
6. The titers of the antibodies on the hemagglutination reactions are not related to the lesions or clinical seriousness.
7. All human components have been tested for anti-HCV, anti-HIV and HBsAg and found negative. However, this does not assure the absence of these and other infectious diseases. All human components should be handled as potentially hazardous..
8. Do not freeze reagents.
9. Sodium azide 0.095% is used as preservative. It can react with lead and cooper plumbing to form highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide buildup. Sodium azide is toxic if ingested.
10. Disposal in accordance with local regulation.
11. Follow the good laboratory practices (GLP) related to storage, handling and material disposal.

WARRANT

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replace if technical support finds evidence that running, handling and storage were not properly followed.

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BIBLIOGRAFIA / BIBLIOGRAPHY / BIBLIOGRAFIA

1. Bing, D.H.; Weyand, J.G.M.; Stavitsky, A.B.: Hemagglutination with aldehyde-fixed erythrocytes for assay of antigens and antibodies. *Proc. Soc. Exp. Bio.*, 124: 1166-1170, 1967.
2. Boné, G. J.; Parent, G.; Stearic acid, an essential growth factor for *Trypanosoma cruzi*. *J. Exp. Microbiol.*, 31: 261-266, 1993.
3. Boyden, S.V.: The adsorption of proteins on erythrocytes treated with tannic acid and subsequent hemagglutination by antiprotein sera. *J. Exp. Med.*, 93: 107-120, 1951.
4. Camargo, M.E.; Hoshino, S.; Siqueira, G.R.V.: Hemagglutination with preserved, sensitized cells: a practical test for routine serologic diagnosis of American trypanosomiasis. *Rev. Inst. Med. Trop. S. Paulo*, 15(2): 81-85, 1973.
5. Camargo, M.E. et al.: Hemagglutination test for Chagas' disease with chromium chloride, formalin-treated erythrocytes, sensitized with *Trypanosoma cruzi* extracts. *Rev. Inst. Med. Trop. S. Paulo*, 13(1): 45-50, 1971.
6. Dias, J.C.P.: Doença de Chagas. In: Cimerman, B. e Cimerman, S.: **Parasitologia humana e seus fundamentos gerais**. Atheneu: 81-111, 1999.
7. Ferreira, A.W. e Avila, S.L.M.: Doença de Chagas. In: Ferreira, A.W.; Avila, S.L.M. (eds). **Diagnóstico laboratorial das principais doenças infecciosas e auto-imunes**. Guanabara Koogan: 114-149, 1996.
8. Hoshino-Shimizu, S. et al.: A control chart method for evaluating hemagglutination reagent used in Chagas' disease diagnosis. **PAHO Bull.**, 20(2): 170-178, 1986.
9. Luquetti, O.A.; Castro, A.M.: Diagnóstico sorológico da doença de Chagas. In: Dias, J.C.P. e Coura, J.R. (eds). **Clínica terapêutica da doença de Chagas. Um manual para o clínico geral**. **Rio de Janeiro:** Fiocruz: 99-114, 1997.
10. Oddo, D. et al.: Acute Chagas' disease (trypanosomiasis americana) in acquired immunodeficiency syndrome. *Hum. Pathol.*, 23: 41-44, 1992.

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