Imuno-HAI TPHA

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ENGLISH

SUMMARY

Syphilis is human infectious disease caused by a spirochete bacterium, *Treponema pallidum*. It is a sexual transmitted disease. Other possible routes of transmission are infected blood transfusion and perinatal (congenital syphilis) from mother to child in utero. Primary syphilis is manifested after an incubation period of 10-90 days. During the initial incubation period, patients are asymptomatic. The sore, called a chancre, is a firm, painless skin ulceration occurring in 85% of patients.

Secondary syphilis is characterized after 4 to 10 weeks after primary infection. The most common symptoms are weight loss, headache, anorexia, myalgia, arthralgia, malaise, fever, lymph nodes and exanthem, occurring in 75 to 100% of patients. In this phase, central nervous system can also be affected. Latent syphilis is after the primary and secondary stages and is characterized by positive serological tests and the absence of clinical findings. This stage can last 1 to 2 years. Tertiary syphilis is manifested in one third of patients without treatment and can present gumma (15%), cardiovascular syphilis (10%) or neuro-syphilis (8-10%). Serological tests for syphilis are classified as non-treponema, for screening, such as VDRL (Veneral Disease Research Laboratory) and RPR (Rapid Plasma Reagin) and treponemal tests, used as confirmatory tests for reactive serum in screening, such as TPHA (*Treponema pallidum* Hemagglutination), FTA-Abs (Fluorescent Treponemal Antibody Absorption) and ELISA (Enzyme-Linked Immunosorbent Assay).

Hemagglutination test for *Treponema pallidum* is sensitive and easy to use.

PRINCIPLE OF THE METHOD

Stabilized fowl erythrocytes sensitized with highly purified *Treponema pallidum* antigens reveal agglutination when react with antibodies against these antigens present in serum or plasma, from patients with suspect of syphilis.

Non-specific antibodies, from treponemal saprophyte that could present false positive results, are absorbed by components in the stabilizer diluent of sensitised erythrocytes.

KIT PRESENTATION

REF 35096-H (96 qualitative tests)

1. Erythrocytes suspension sensitised with components of *Treponema pallidum* (1x5ml)
2. Diluent solution (1x40ml)
3. Positive control serum (1ml)
4. Negative control serum (1ml)
5. Disposable "V"-well microtitration plates (1x96)
6. Instructions for use.

REF 35192-H (192 qualitative tests)

1. Erythrocytes suspension sensitised with components of *Treponema pallidum* (2x5ml)
2. Diluent solution (2 x 40ml)
3. Positive control serum (1ml)
4. Negative control serum (1ml)
5. Disposable "V"-well microtitration plates (2x96)
6. Instructions for use.

REAGENT STABILITY AND STORAGE

• **ERYTHROCYTES SUSPENSION (1):** Ready for use. Stable if stored at 2-8°C up to expiration date. **Do not freeze.** Homogenize before using it, always in circular movements and for vial inversion but never shake, since this movement can induce to haemolysis and/or agglutination. Allow the suspension to reach room temperature prior to use.

• **DILUENT SOLUTION (2):** Ready for use. Stable if stored at 2-8°C. Allow the components to reach room temperature prior to use. It contains Sodium azide 0.095%.

• **POSITIVE CONTROL SERUM (3):** Ready for use. Stable if stored at 2 to 8°C.

• **NEGATIVE SERUM CONTROL (4):** Ready for use. Stable if stored at 2 to 8°C.

The kit presents good performance after being used for the first time. It is stable up to the expiration date if store at 2-8°C.

SPECIMEN COLLECTION AND STORAGE

Fresh serum or plasma from patients or stored at -20°C (4 to 6 weeks) should be used. The specimen should reach to room temperature prior to use. Do not use haemolysed, contaminated or lipaemic specimen.

MATERIAL REQUIRED BUT NOT PROVIDED

- Test tubes for dilution and titration
- Serological pipettes
- Shelf for tubes and rack for tips
- Container for discarded material

PROCEDURE

A. Qualitative test

1. Lay the plate onto a humid piece of cloth to avoid electrostatics.
2. Use 1 well of the plate for specimen, including positive and negative controls.

- Dilute in the test tube: 390 μ l of diluent (2) + 10 μ l of specimen (dilution at 1:40).
- Add 50 μ l of positive control (3), negative control (4) in the A1 and A2 wells, for example. (ATTENTION: DO NOT DILUTE THE CONTROLS. THEY ARE READY FOR USE). Add 50 μ l of the dilution at 1:40 of each specimen in the respective wells of the plate.
- Add 50 μ l of the homogeneous suspension of erythrocytes (1) on each well, which correspond to one dilution of the specimen (unknown serum) at 1:80.
- Shake the plate by mechanical vibration (plate agitator) or beat the fingers on the edges of the plate for 3 to 4 minutes.
- Let stand for 1 hour at room temperature in a place free of vibration (IMPORTANT).
- Interpret the results.

RESULTS

Negative: Red cells gather in the bottom of the well and form a button.

Positive: Red cells gather in the bottom of the well and form a rug, sometimes it shows irregular edges.

INTERPRETATION

The agglutination of test cells indicates the presence of specific antibody to anti- **Treponema pallidum**. The absence of agglutination indicates the absence of specific antibody to anti-**Treponema pallidum** or antibody is below the limit of detection of the system.

NOTE: Hemagglutination test associated with other methods such as immunofluorescence and/or ELISA are recommended as confirmatory tests.

B. Semi-Quantitative test

- Dilute the specimen at 1:40 as per step 3 of the qualitative test.
- Add 100 μ l of the specimen diluted at 1:40 in the well A1 of the plate.
- Add 50 μ l of the diluent solution in the wells B1 to H1 of the plate.
- Transfer 50 μ l from the well A1 to the well B1 by using the automatic pipette and mix gently. Transfer 50 μ l from well B1 to well C1 and so on up to well H1. Discard 50 μ l of the last well.
- Add 50 μ l of the suspension of erythrocytes in all wells which contain specimens. The dilutions vary from 1:80 (A1) to 1:10240 (H1).
- Shake the plate by mechanical vibration (plate agitator) or beat the fingers on the edges of the plate during 3 to 4 minutes.
- Let stand for 1 hour at room temperature in a place free of vibration (IMPORTANT).
- Interpret the results.

TITER OF THE SPECIMEN: The highest dilution which shows a positive reaction. The end point is when the erythrocytes cover 50% of the bottom of the well.

PERFORMANCE OF THE TEST

Sensitivity: 100% of sensitivity - 95 in-house positive specimens from the quality control were compared. All results were satisfactory. False negative results were not found.

Specificity: 99,7% of specificity - 275 in-house negative specimens from the quality control were compared. It was found one false positive result.

Precision: Intra-assay: The precision was verified using 10 replicas of 10 specimens. The negative and positive results were correctly identified at 100% of the times.

Inter-assay: The precision was verified using 10 replicas of 10 specimens against independent assay with 3 different lots during a period of 3 days. The negative and positive results were correctly identified at 100% of the times.

PRECAUTIONS AND WARNINGS

- Do not use reagents from different lots.
- The plates should be free of vibration since it can affect the interpretation of the results
- The vial of the erythrocytes should be stored in vertical position to avoid irreversible auto-agglutination.
- USE EXCLUSIVELY THE KIT PLATE. It should not be reused.
- Variations of 1 or 2 titers are acceptable for the same serum on hemagglutination reactions.
- No serological test through hemagglutination method can distinguish antibodies due to **Treponema pallidum** infection because of other pathogenic treponemal infections such as **T. pertenuae** and **T. Carateum**.
- A confirmatory FTA Abs test is also recommended, kit **Imuno ConSifilis (FTA-Abs)** from **WAMA**.

8. Test can be negative in early active syphilis or later latent syphilis. To complete and aid the syphilis diagnostic, it is recommended the use of **VDRL** or **RPR**. **WAMA** markets both tests, **VDRL** and **RPR** (Reagins).

9. All human components have been tested for anti-HCV, anti-HIV and HBsAg and found negative. However, this does not assure the absence of these and other infectious diseases. All human components should be handled as potentially hazardous.

10. Do not freeze the reagents.

11. Sodium azide 0.095% is used as preservative. It can react with lead and copper plumbing to form highly explosive metal azides. Upon

disposal, flush with large volumes of water to prevent azide buildup. Sodium azide is toxic if ingested.

12. Disposal in accordance with local regulations.

13. Follow the good laboratory practices (GLP) related to storage, handling and material disposal.

WARRANT

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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