Imuno-HAI TOXOPLASMOSE

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ENGLISH

SUMMARY

Toxoplasmosis is caused by a protozoan called *Toxoplasma gondii*. As this protozoan can parasite all nucleated cells, the infection has a wide range of impact, such as chorioretinitis, myocarditis, pericarditis, pneumonia, myositis, meningoencephalitis, hepatitis; these infections can be isolated or combined, often diluting, particularly in patients with immunodeficiency or undergoing immunosuppressant therapy. The most serious problem of *T. gondii* infection is the congenital transmission during pregnancy. The severe compromise of the fetus may lead to the death of the fetus, death of the newborn or congenital disease, expressed by the following symptoms: hydrocephalus or microcephalus, Chorioretinitis, and brain calcification. The transplacental congenital transmission of the toxoplasmosis depends on the existence of the acute infection during the gestational phase; the earlier the gestational age, the greater the risk, being the highest risk during the first trimester.

Usually the beginning of the *T. gondii* infection is insidious, characterized by fever, adynamia, adenopathy, especially the cervical one; however, it frequently manifests in subclinical form, presenting in unpecific and light symptoms.

The toxoplasmosis diagnosis should never be made solely on the findings of clinical data. There is always the need to confirm it through laboratory tests and correct interpretation of the results.

The most important laboratory methods for the diagnosis of toxoplasmosis are: indirect immunofluorescence, passive hemaglutination and ELISA.

PRINCIPLE OF THE METHOD

Fowl erythrocytes sensitised and stabilized with highly purified *Toxoplasma gondii* antigens show agglutination when react with antibodies against these antigens present in patient's serum.

KIT PRESENTATION

REF 36096-H (96 qualitative determinations)

1. Suspension of sensitised erythrocytes with *Toxoplasma gondii* (2.4ml)

2. Diluent solution (1x 40ml)

3. 2-Mercaptoethanol (1 x 0.5ml)

4. Positive serum control (1 x 1ml)

5. Negative serum control (1 x 1ml)

6. Disposable “V”-well microtitration plate (1x96 wells)

7. Instructions for use

REF 36192-H (192 qualitative determinations)

1. Suspension of sensitised erythrocytes with *Toxoplasma gondii* (1 x 4.8ml)

2. Diluent solution (1x60ml)

3. 2-Mercaptoethanol (1 x 0.5ml)

4. Positive serum control (1 x 1ml)

5. Negative serum control (1 x 1ml)

6. Disposable “V”-well microtitration plate (2x96 wells)

7. Instructions for use

REF 36380-H (380 qualitative determinations)

1. Suspension of sensitised erythrocytes with *Toxoplasma gondii* (2x4.8ml)

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7. Instructions for use

REAGENT PREPARATION AND STABILITY

• **SUSPENSION OF SENSITISED ERYTHROCYTES (1):** Ready for use. Keep at 2-8°C. Do not freeze.

Homogenize the erythrocytes prior to use by inverting and round movements, never shake the vial since it can cause haemolysis and/or auto-agglutination.

• **DILUENT SOLUTION (2):** Ready for use. Keep at 2-30°C.

• **2-MERCAPTOETHANOL (3):** Protect from light. Dilute with the diluent solution (2). Ex: 10ml of diluent solution for 70µl of 2-mercaptoethanol (3). This mixture is intended to dilute the serum in order to determine a recent infection (IgM) or the presence of non-specific antibodies.

The solution can be diluted in ambar vial and stored at 2-8°C together with the kit. After preparing the diluent solution with 2-ME is stable for 30 days.

• **POSITIVE SERUM CONTROL (4):** Ready for use. Keep at 2-8°C.

• **NEGATIVE SERUM CONTROL (5):** Ready for use. Keep at 2-8°C.

The kit presents good performance after used for the first time. It is stable up to the expiration date if store at 2-8°C.

SPECIMEN COLLECTION

Use fresh serum or stored at -20°C (maximum of 4 to 6 weeks) from fasting patients. The specimens must be at room temperature before using and do not need to be treated. Do not use hemolized, contaminated or lipemic serum.

Freezing and thawing cycles should be avoided. Old serum tend to become gelatinous in contact with 2-mercaptoethanol.

MATERIAL REQUIRED BUT NOT PROVIDED

- Test tubes
- Serological pipettes
- Tube support and rack
- Disposable material container

PROCEDURE

a. QUALITATIVE TEST (“Screening”)

Objective: Intended for screening and ruling out of non-reactive serum.

1. Lay the plate onto an humid cloth to avoid electrostatic.

2. Use 1 well for each specimen, including positive and negative controls.

ATTENTION: Do not dilute the controls. They are ready for use.

3. Dilute the serum 1/32 with diluent solution (2). It's recommended to dilute 10µl of

specimen + 310µl of diluent solution in a test tube.

4. Pipette 25µl of positive and negative controls and serum diluted at 1/32 of each sample on respective wells. Suggestion: A1 positive control, A2 negative control and A3, A4, A5...sera to be tested.

5. Add 25µl of suspension of erythrocytes (1) into each well.

6. Shake the plate by mechanical vibration (plate agitator) or beating with the fingers on the edges of the plate for 3 to 4 minutes.

7. Rest the plate for 1 to 2 hours in a place free of vibrations at room temperature (IMPORTANT).

8. Read the results. (See semi quantitative procedure).

INTERPRETATION

Serum that has a negative reaction are reported as NON-REACTIVE; they do not have current or previous toxoplasmosis infection. When one observes a negative reaction, but suspects toxoplasmosis, the test should be repeated after a few days since it may fall within the pre-serological period. All sera that gives positive reactions should be quantitatively tested to determine the stage of the disease

B. SEMI-QUANTITATIVE TEST

1. Dilute the serum 1/32 according to step 3 of qualitative test.

2. Pipette 25µl of diluent solution (2) from the second well of the plate up to dilution that is intended to evaluate. Ex: If the dilution of the serum is up to 1/512, pipette 25µl of diluent in the wells A2, A3, A4 and A5. Do not pipette in the well A1.

3. Transfer 25µl of the serum diluted at 1/32 in the wells A1 and A2 of the plate. Homogenize thoroughly the serum with diluent in the well A2 (dilution 1/64) and transfer 25µl in the well A3 (dilution 1/128) and so on. Discard 25µl in the end.

4. Pipette 25µl of erythrocytes homogeneous suspension (1) in all wells.

5. Shake the plate by mechanical vibration (plate vibrator) or by beating with the fingers on the edges of the plate for 3 to 4 minutes.

6. Rest the plate for 1 to 2 hours in a place free of vibrations at room temperature (IMPORANT).

7. Read the results.

READING

Negative: Erythrocytes form a compact button at the bottom of the well.

Positive: Erythrocytes form an even layer at the bottom of the well, like a rug, often presenting irregular edges.

TITER OF THE SPECIMEN: The highest dilution showing positive reaction. The endpoint is when the erythrocytes cover 50% of the bottom of the well.

DETERMINATION OF RECENT INFECTION (IgM) OR UNSPECIFIC ANTIBODIES

The use of 2- mercaptoethanol, on the hemagglutination tests, offers the possibility to differentiate the recent infection or the presence of non-specific antibodies.

PROCEDURE

1. Perform the test with 2-Mercaptoethanol in parallel to Quantitative Test.

2. Serum treatment: Pipette 10µl of the activated or inactivated serum + 310µl of diluted 2-Mercaptoethanol (1/32 solution) in a test tube. Incubate at 37°C for 1 hour.

INTERPRETATION

- Negative: Patients non-immune or in pre-serological stage. When in doubt, repeat the test after 1 or 2 weeks.

- Titers of serum that have been treated and non-treated (with 2-mercaptoethanol) that are similar (less than 2 dilutions): Probably previous toxoplasmosis infection.

ATTENTION: This result does not exclude recent infection (IgM) and therefore, immunofluorescence and ELISA for IgM detection are recommended.

- Difference of 2 or more titers for the treated serum versus the non-treated (with 2-mercaptoethanol): Probably recent infection (IgM)

ATTENTION: If titers are low and treated serum is negative, to rule out the presence of unspecific antibodies, repeat the test after 1 to 2 weeks and observe if the titer is higher.

NOTE: The hemagglutination test associated with more sensitive methods such as immunofluorescence and/or ELISA enables one to establish serological profiles that best characterizes the disease stage.

PERFORMANCE CHARACTERISTICS

In a study on 74 true reactive specimens performed by indirect immunofluorescence and ELISA methods, the **Imuno-HAI TOXOPLASMOSE** from **WAMA** Diagnóstica showed no false negative results which represents a sensitivity of 100%.

Another study on 268 true negative specimens performed by indirect immunofluorescence and ELISA methods, the **Imuno-HAI TOXOPLASMOSE** from **WAMA** Diagnóstica showed 2 false positive results which represents a specificity of 99.3%.

PRECAUTIONS AND WARNINGS

1. In vitro diagnostic.

2. USE ONLY THE PLATE THAT COMES WITH THE KIT PLATE. It should not be reused.

3. Do not use reagents from different lots.

4. When reading the plate, the plate should not vibrate since it can affect the reading results.

5. Erythrocytes vial should be let stand in vertical position to avoid irreversible auto-agglutination.

6. It is acceptable variations of 1 or 2 titers for the same serum on hemagglutination reactions.

7. All human components used in the preparation of the controls have been tested (and found negative) for anti-HIV, HBsAg, and anti-HCV. However, since no diagnostic method can give full assurance of the absence of these and other infectious agents, it's recommended to treat the human serum controls as potentially infective materials

8. Do not freeze the reagents.

9. Sodium azide 0.095% is used as preservative. Upon disposal, flush with large volumes of water to prevent azide buildup in the plumbing, which can become explosive.

10. Disposal in accordance with local regulations.

11. Follow good laboratory practices (GLP) related to storage, handling and material disposal.

WARRANTY

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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