MS 10310030077
10310030088*

Imuno-Látex

FR

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ENGLISH

SUMMARY

The diagnosis of Rheumatoid arthritis is based on clinical assay but the laboratory and radiological tests are useful to confirm clinical diagnosis and to evaluate the seriousness and the phase of the disease. Rheumatoid Factor (RF) in serum is one of the main clinical markers for rheumatoid arthritis. RF is a term used to describe a group of antibodies (IgM, IgG, IgA and IgE) directed

against the antigenic site on the Fc portion of the body's own IgG antibodies.

Many tests have been suggested for the detection of RF such as the ones that use human erythrocytes coated with human or animal IgG. Other methods, more sensitive, have also been described such as latex particles coated with human IgG.

Imuno Latex FR is a very sensitive test using a suspension of polystyrene latex particles coated with purified and stabilized human IgG and suspended in glycine buffer pH 8.2.

PRINCIPLE OF THE METHOD

A suspension of latex is mixed on an test card with serum containing elevated rheumatoid factor levels, a clear agglutination is seen in 2 minutes.

KIT PRESENTATION

REF 2860-L (60 determinations)

1. Suspension of latex (1.5ml)
2. Positive control serum (0.5ml)
3. Negative control serum (1ml)
4. Stirrers (60)
5. Test card (2)
6. Instructions for use

REF 28100-L (100 determinations)

1. Suspension of latex (2.5ml)
2. Positive control serum (0.5ml)
3. Negative control serum (1ml)
4. Stirrers (100)
5. Test card (2)
6. Instructions for use

REF 2800-L (100 determinations)

1. Suspension of latex (2.5ml)
 2. Instructions for use
- MATERIAL REQUIRED BUT NOT PROVIDED**
- Test tube
 - Serological pipettes
 - Shelf for tubes and rack
 - Disposable material container
 - Saline 0.9%

REAGENT STABILITY AND STORAGE

• **SUSPENSION OF LATEX (1):** Stable if stored at 2-8°C up to expiration date. Allow the suspension to reach to room temperature and mix gently prior to use. Do not freeze. It contains Sodium azide 0.1%.

• **POSITIVE SERUM CONTROL (2):** Ready for use. Stable if stored at 2-8°C up to expiration date. Allow the control to reach to room temperature prior to use. RF concentration: > or = 8IU/ml. It contains Sodium azide 0.1%.

• **NEGATIVE SERUM CONTROL (3):** Ready for use. Stable if stored at 2-8°C up to expiration date. Allow the control to reach to room temperature prior to use. It contains Sodium azide 0.1%.

The kit presents good performance after being used for the first time. It is stable up to the expiration date if store at 2-8°C.

SPECIMEN COLLECTION AND STORAGE

Do not use Haemolysed, lipaemic or contaminated serum for testing. The specimen may be stored at 2-8°C for up to 48 hours prior testing. If longer storage is required, store it in freezer at -20°C for up to 6 weeks. Undiluted samples should be used. Do not use plasma since fibrinogen can form non-specific agglutination

PROCEDURE

A. Qualitative test

1. Transfer 25µl of patient's serum to the test circle of the test card.
2. Shake gently the suspension of latex (1), then using the pipette, add 25µl of the suspension to the same test circle.
3. Using the stirrers, mix the serum and the latex and spread them over the entire circle.
4. Rotate the test card for 2 minutes and observe the agglutination while holding the test card under a good source of light.

ATTENTION: Controls should be tested with each test run in

order to check the procedure and reagents stability. Prozone effect can be present in concentrations higher than 1000IU/ml. In case of higher values, the specimen should be diluted.

RESULTS

Positive: It is indicated by the obvious agglutination pattern of the latex. RF serum concentration of 8IU/ml or more.
Negative: It is indicated by no change in the latex suspension on the test card. RF concentration below 8IU/ml.

ATTENTION: The sensitivity of the test is 8IU/ml. Therefore, negative serum is below 8IU/ml. The sensitivity has been adjusted according to the International Standard for the arthritic serum from WHO.

B. Semi-Quantitative test

1. Using isotonic saline (NaCl 0.9%) prepare serial dilutions of the patient's serum at 1:2, 1:4, 1:8, 1:16 and so on.

2. Transfer 25µl of each serum dilution to the test circle of the test card.

3. Mix suspension of latex (1), then using a pipette, add 25µl of the suspension to the test circle where the serum dilution is present.

4. Using the stirrer, mix the serum with the latex and spread them over the entire circle (Use a separate stirrer for each dilution).

5. Rotate the test card for 2 minutes and observe the agglutination while holding the slide under a good source of light.

The highest dilution of the serum which shows an obvious agglutination pattern is considered the titer end point.

The serum RF concentration is calculated approximately by multiplying the dilution factor (i.e. 2,4,8 or 16) by the detection limit, i.e. 8, to give the number of IU/ml concentration.

EXAMPLE: If the agglutination titer appears at 1:8, the approximate serum RF concentration is: $8 \times 8 = 64$ IU/ml.

INTERPRETATION

RF positive is present between 10 and 80% of patients with rheumatoid arthritis. However, they are not specific of arthritis once patients suffering from syphilis, LES, infectious mononucleosis, hepatitis, hypergammaglobulinemia etc, can present positive results, but with low titers. Less than 5% of healthy patients can present positive RF. On the other hand, some patients with rheumatoid arthritis can present RF negative. RF positive is usually present after 6 months of the onset of the disease.

PRECAUTIONS AND WARNINGS

1. Store the reagents at 2-8°C. Do not freeze.

2. Read the results after 2 minutes can cause false positive results.

3. At the end of the test, wash the test card with distilled water. If prefer to wash later, use water with a neutral detergent and rinse out slightly with distilled or deionized water. Dry the test card prior to reuse.

4. Use a separate stirrer for each specimen.

5. The reagents (latex and controls) contain 0.1% sodium azide as a preservative which may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. On disposal, flush with large quantities of water.

6. The source material used in the preparation of the controls has been tested and confirmed negative for HCV, HIV and HBsAg. Because no test can offer complete assurance that products derived from human source will not transmit infectious agents it is recommended that these controls be treated as potential biohazards materials.

7. False positive reactions which are found by using latex methods, due to the presence of complement (C1q), can be avoided by heat inactivation or specimen dilution.

8. Disposal in accordance with local regulations.

9. Follow the good laboratory practices (GLP) related to storage, handling and material disposal.

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