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## ENGLISH

### SUMMARY

Infectious Mononucleosis (IM) is an acute disease transmitted by Epstein-Barr virus. Mostly benign in children and teenagers, this disease presents posterior cervical adenopathy fever and throat inflammation.

IM also presents an increase of monocytes and lymphocytes with more than 10% of different forms. Serologically, IM presents agglutinin of sheep cells, haemolysin of bovine cells and the presence of antibodies against Epstein-Barr virus. In 1932, Paul and Bunnell were the first to report that antibodies from IgG class, which agglutinated the erythrocytes of sheep and horse, were present in the patient's serum. Later on, Davidson employed a differential absorption procedure using kidney extractor of cavy or horse.

Imuno-Latex MNI from WAMA Diagnostica presents a suspension of latex particles coated with antigenic extract of stroma of erythrocytes from bovine.

The specificity of the test is obtained through the use of purified antigens avoiding the absorption of interfering antibodies. The test is simple and easy to read.

### PRINCIPLE OF THE METHOD

A macroscopic agglutination occurs when the serum or plasma containing heterophile antibodies associated with mononucleosis are mixed with latex reagent.

### KIT PRESENTATION

REF 21160-L (60 determinations)

1. Suspension of latex coated with purified antigens from bovine erythrocytes (1,5ml)
2. Positive serum control (0,5ml)
3. Negative serum control (1ml)
4. Stirrers (60)
5. Test card (2)
6. Instructions for use

REF 211100-L (100 determinations)

1. Suspension of latex coated with purified antigens from bovine erythrocytes (2,5ml)
2. Positive serum control (0,5ml)
3. Negative serum control (1ml)
4. Stirrers (100)
5. Test card (2)
6. Instructions for use

### MATERIAL REQUIRED BUT NOT PROVIDED

- Test tube
- Serological pipettes
- Tube shelf and rack for tips
- Disposable material container
- Saline 0.9%

### STABILITY AND STORAGE

• **SUSPENSION OF LATEX (1):** Stable if stored at 2-8°C up to expiration date. Allow the latex to reach room temperature and mix it prior to use. Do not freeze. It contains sodium azide 0.1%

• **POSITIVE SERUM CONTROL (2):** Ready for use. Stable if stored at 2-8°C up to expiration date. Allow the control to reach room temperature prior to use. It contains sodium azide 0.1%.

• **NEGATIVE SERUM CONTROL (3):** Ready for use. Stable if stored at 2-8°C up to expiration date. Allow the control to reach room temperature prior to use. It contains sodium azide 0.1%.

**Obs.:** The kit presents good performance after being used for the first time. It's stable up to the expiration date if store at 2-8°C.

#### SPECIMEN COLLECTION AND STORAGE

Serum or plasma from blood with heparin or EDTA. Haemolysis, lipaemic or contaminated serum must be discarded. The specimen may be stored at 2-8°C for a week prior testing. If longer storage is required, store in freezer at -20°C for up to 6 weeks.

Undiluted serum or plasma is recommended.

#### PROCEDURE

1. Pipette 25µl of patient's serum in the area of test card.
2. Gently shake suspension of latex (1) and pipette 25µl in the same area of test card.
3. Using the stirrers, mix the serum dilution and the latex and spread them over the entire circle.
4. Rotate the test card for 2 minutes and observe the agglutination while holding the slide under a good source of light.

#### ATTENTION:

1) Positive and negative controls should be included at regular intervals in order to check the procedure and reagents.

2) If patient presents mononucleosis but the test result is negative, repeat the test after diluting the specimen at 1:10 with saline (0.9% Na Cl) since the Prozone effect can be present.

#### READING

Positive: Clear latex agglutination is seen.

Negative: No latex agglutination is seen.

#### PRECAUTIONS AND WARNINGS

1. Store reagents at 2-8°C. Do not freeze.
2. At the end of the test, wash the test card with distilled water. If prefer to wash later, use water with a neutral detergent and rinse out slightly with distilled or deionized water. Dry prior to use.
3. Use separate stirrer for each sample
4. Imuno Látex MNI (latex and controls) contains 0.1% sodium azide as a preservative which may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. On disposal, flush with large quantities of water.
5. The source material used in the preparation of the controls has been tested and confirmed negative for HCV, HIV and HBsAg. Because no test can offer complete assurance that products derived from human source will not transmit infectious agents it is recommended that these controls be treated as potential biohazards materials.
6. Disposal in accordance with local regulations.
7. Follow the good laboratory practices (GLP) related to storage, handling and material disposal.

#### WARRANTY

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by

Wama's technical support. The warranty will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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