MS 10310030072
10310030090*

Imuno-Látex PCR

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ENGLISH

SUMMARY

C-Reactive Protein together with Erythrocyte Sedimentation Rate (ESR) and mucroprotein are known as "acute infections". They characterize non-specific infections, helping the diagnosis and controlling the development of the infection.

C-Reactive Protein (CRP) in patient's serum has been found in association with acute infections, necrotic conditions such as bacteria and virus infections, acute rheumatic fever, myocardial infarction, tuberculosis, rheumatoid arthritis and malign neoplasia. The increase in the level of CRP is a non-specific answer to the inflammatory process.

Imuno-Látex CPR contains latex particles coated with monoclonal anti-CRP antibody, stabilized and suspended in glycine-buffer, pH 8.2.

PRINCIPLE OF THE METHOD

Suspension of latex is mixed with serum containing CRP on a circle of test card and a clear agglutination is seen if the CRP concentration is from 6 to 400mg/L. Agglutination is not seen if CRP concentration is below 6mg/L.

KIT PRESENTATION

REF 2960-L (60 determinations)

1. Suspension of latex coated with monoclonal anti-PCR antibody (1.5ml).
2. Positive serum control (0.5ml)
3. Negative serum control (1.0ml)
4. Stirrers (60)
5. Test cards (2)
6. Instructions for use

REF 29100-L (100 determinations)

1. Suspension of latex coated with monoclonal anti-PCR antibody (2.5 ml).
2. Positive serum control (0.5ml)
3. Negative serum control (1.0ml)
4. Stirrers (100)
5. Test cards (2)
6. Instructions for use

REF 2900-L (100 determinations)

1. Suspension of latex coated with monoclonal anti-PCR antibody (2.5 ml).
 2. Instructions for use
- MATERIAL REQUIRED BUT NOT PROVIDED**
- Test tube 13 X 75 mm
 - Serological pipettes
 - Test tube shelf and racks for tips
 - Disposable material container
 - Saline 0.9%

REAGENT STABILITY AND STORAGE

• **SUSPENSION OF LATEX (1):** Stable if stored at 2-8°C up to expiration date. Allow the suspension to reach room temperature and mix it prior to use. Do not freeze. It contains Sodium azide 0.1%.

• **POSITIVE SERUM CONTROL (2):** Ready for use. Stable if stored at 2-8°C up to expiration date. Allow the serum to reach room temperature prior to use. PCR concentration: > or = 10mg/L. It contains Sodium azide 0.1%.

• **NEGATIVE SERUM CONTROL (3):** Ready for use. Stable if stored at 2-8°C up to expiration date. Allow the serum to reach room temperature prior to use. It contains Sodium azide 0.1%. The kit present good performance after being used for the first time. It is stable up to the expiration date if stored at 2 - 8°C.

SPECIMEN COLLECTION AND STORAGE

Haematic, lipaemic or contaminated serum must be discarded. The specimen may be stored at 2-8°C for up to 48 hours prior testing. If longer storage is required, store in freezer at -20°C for up to 6 weeks. Undiluted specimen should be used. Do not use plasma since the fibrinogen can form non-specific agglutination.

PROCEDURE

A. Qualitative test

1. Pipette 25µl of patient's serum in the area of test card.
2. Gently shake the suspension of latex (1) and pipette 25µl in the same area of test card.
3. Using the stirrers, mix the serum and the latex and spread them over the entire circle.
4. Rotate the test card for 2 minutes and observe the agglutination while holding the slide under a good source of light.

ATTENTION: Positive and negative controls should be included at regular intervals in order to check the procedure and reagents.

Prozone effect can be present in concentrations higher than 400mg/L. If it is the case, the sample should be diluted.

READING

Positive:

Clear agglutination is seen. Concentration equal or > 6mg/L.

Negative:

No agglutination is seen. Concentration < 6mg/L.

ATTENTION: Sensitivity is 6mg/L. Therefore, negative serum is lower than 6mg/L.

B. Semi-Quantitative test

1. Using saline (Na Cl 0.9%) prepare serial dilutions of the patient's serum at 1:2, 1:4, 1:8, 1:16 and so on.
2. Pipette 25µl of each dilution in each area of test card.
3. Mix the suspension of latex (1) and pipette 25µl in each area where the sample dilution is present.
4. Using the stirrer, mix the serum dilution with the latex and spread over the entire circle (Use a separate stirrer for each dilution).
5. Rotate the test card for 2 minutes and observe the agglutination while holding the slide under a good source of light. (See Reading).

The highest serum dilution showing a definite agglutination pattern is considered the titer end point. PCR concentration is calculated as follow:

Concentration (mg/L) = 6 X D, where 6 is the sensitivity of the test and D is the highest dilution showing agglutination.

EXAMPLE: If the titer is 1:8, the concentration of PCR in the specimen is: 6x8 = 48mg/L.

INTERPRETATION

Most normal, healthy humans have serum CRP concentrations equal or > 6mg/L. At the onset of disease, plasmatic PCR concentrations increase rapidly when compared with erythrocyte sedimentation rate and it is not seen after the patient's recovering.

CRP detection can be valuable as an aid in monitoring the effectiveness of therapy and inflammatory conditions particularly in the management of patients with acute rheumatic fever. When CRP levels, under treatment, reappears, it is considered a worst on the treatment process. In these situations, it is recommended a monitoring at regular time.

PCR post-surgery monitoring is an aid as an indicator of septic or inflammatory state. PCR levels can be increased at the first three months of pregnancy and remain up high to parturition. The levels also are increased in women using oral contraceptives.

PERFORMANCE CHARACTERISTICS

Imuno-Rápido PCR presented 100% of sensibility and 100% of specificity. 50 specimen divided in 13 positive and 37 negative were tested and confirmed against kits from the market. It was not observed neither false positive nor false positive results.

PRECAUTIONS AND WARNINGS

1. Store reagents at 2-8°C. Do not freeze.
2. Read the results after 2 minutes can cause false positive results.
3. At the end of the test, wash the test card with distilled water. If prefer to wash later, use water with a neutral detergent and rinse out slightly with distilled or deionized water.
4. Use separate stirrer for each sample.
5. Imuno-Látex PCR (latex and controls) contains 0.1% sodium azide as a preservative which may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. On disposal, flush with large quantities of water.

6. The source material used in the preparation of the controls has been tested and confirmed negative for HCV, HIV and HBsAg. Because no test can offer complete assurance that products derived from human source will not transmit infectious agents it is recommended that these controls be treated as potential biohazards materials.

7. Disposal in accordance with local regulations.

8. Follow the good laboratory practices (GLP) related to storage, handling and material disposal.

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