RPR Sífilis

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ENGLISH

SUMMARY

Syphilis is a sexually transmitted disease (STD) caused by the spirochete, *Treponema pallidum*.

Other transmission ways are infected blood transfusion, which is rare today due to serological screening, and perinatal (congenital syphilis) that is transmitted from mother to child in utero. Primary syphilis is manifested after an incubation period of 10-90 days. During the initial incubation period, patients are asymptomatic. The sore, called chancre, is a firm, painless skin ulceration occurring in 85% of patients.

Secondary syphilis is characterized after 4 to 10 weeks after primary infection. The most common symptoms are weight loss, headache, anorexia, myalgia, arthralgia, malaise, low fever,

limph nodes and exanthem, occurring in 75 to 100% of patients. On this phase, central nervous system can also be affected. Latent syphilis occurs after the primary and secondary stages and is characterized by positive serological tests and the absence of clinical findings. This stage can last 1 to 2 years. Tertiary syphilis is manifested in one third of patients without treatment and can present gumma (15%), cardiovascular syphilis (10%) or neuro-syphilis (8-10%)

Serological tests for syphilis are classified as non-treponemal such as **VDRL (Venereal Disease Research Laboratory)** and **RPR (Rapid Plasma Reagin)** and treponemal tests such as **TPHA (Treponema pallidum Hemagglutination)**, **FTA-Abs (Fluorescent Treponemal Antibody Absorption)**.

RPR (Rapid Plasma Reagin) is a flocculation, non-treponemal test used to detect syphilis through research of antibodies (reagin) in serum or plasma. It is a modified form of the classic VDRL.

RPR contains coal particles enabling a visual reading of results and avoiding microscopy devices, which is necessary for VDRL.

The test can be performed in serum or plasma, activated or inactivated. RPR kit is provided with antigenic suspension, test cards, stirrers, positive and negative controls.

PRINCIPLE OF THE METHOD

When RPR antigenic suspension is mixed with serum or plasma which contain antibodies (reagin), the particles of antigens flocculate and the result of the reaction is seen under a good source of light. The absence of flocculation indicates a negative result.

KIT PRESENTATION

REF 55100-F (100 determinations)

1. Antigenic suspension (2.5ml)
2. Positive serum control (1ml)
3. Negative serum control (1ml)
4. Stirrers (100)
5. Test card (2)
6. Package insert

REF 55200-F (200 determinations)

1. Antigenic suspension (5ml)
2. Positive serum control (1ml)
3. Negative serum control (1ml)
4. Stirrers (200)
5. Test card (4)
6. Package insert

MATERIAL REQUIRED BUT NOT PROVIDED

- Test tube
- Serological pipettes
- Test tube rack and tip rack
- Disposable material container
- Saline 0.9%

SPECIMEN COLLECTION AND HANDLING

Only serum or plasma should be used with this test, activated or inactivated. Specimen free of haemolysis, lipaemy and bacterial contamination is recommended. Specimen can be stored at -20°C for 4 to 6 weeks, if needed.

STABILITY AND STORAGE

• **ANTIGENIC SUSPENSION (1):** Bring the suspension to room temperature and mix gently prior to use. Stable if stored at 2-8°C up to expiration date. **Do not freeze.**

• **POSITIVE SERUM CONTROL:** Ready for use. Stable if stored at 2-8°C up to expiration date. Contains Sodium azide 0.095%. Allow the control to equilibrate to room temperature prior to use.

• **NEGATIVE SERUM CONTROL:** Ready for use. Stable if stored at 2-8°C up to expiration date. Contains Sodium azide 0.095%. Allow the control to equilibrate to room temperature prior to use.

The kit presents good performance after being used for the first time. It is stable up to the expiration date if stored at 2-8°C.

PROCEDURE

a. Qualitative test

Objective: Screening and removal of negative specimens.

1. Pipette 25µl of patient's serum or plasma on test card circle.
2. Using a stirrer, spread the specimen over the entire circle.
3. Pipette 25µl of antigenic suspension on the specimen to be tested.
4. Using a mechanical rotator set at 100rpm, shake test card for 8 minutes.
5. Immediately after 8 minutes, examine, under a good source of light, the presence or not of aggregates.

ATTENTION: Positive and negative controls should be run at the same time to keep testing conditions the same and to monitor the performance of the reagent.

READING

Negative: ABSENCE of aggregate. Greyish homogeneous aspect.

Weak Positive: PRESENCE of small dispersed aggregates.

Positive: PRESENCE of medium and large aggregates.

b. Semi-Quantitative test

1. Dilute the specimen in a saline solution at 1:2, 1:4, 1:8, 1:16 and more, if needed.
2. Pipette 25µl of each dilution on each test card circle.
3. Using a stirrer, spread the dilutions over the entire circle.
4. Pipette 25µl of antigenic suspension on each area. Use a separate stirrer for each dilution.
5. Using a mechanical rotator set at 100rpm, shake the test card for 8 minutes.
6. Immediately after 8 minutes, examine, under a good source of light, the presence of aggregates in each dilution.

TITER OF SAMPLE

It will be the last dilution titer in which the presence of aggregates can be seen.

INTERPRETATION

False positive reactions with RPR can occur in some conditions such as: immunizations, infections, pregnancy, malaria, autoimmune diseases (SLE etc), malign diseases etc.

In case of RPR is positive, specimens should be investigated by a confirmatory and specific treponemal test. **WAMA Diagnostica** manufactures and markets **FTA-Abs** (Immunofluorescence) and **TPHA** (Hemagglutination) kits.

TEST PERFORMANCE

RPR-SÍFILIS from **WAMA Diagnóstica** presented a sensitivity of 100%. 139 known positive specimens were used. In 203 known negative specimens, any false positive result was found that confirmed to RPR-SÍFILIS from **WAMA Diagnóstica** an specificity of 100%.

PRECAUTIONS AND WARNINGS

1. Store the reagents at 2-8°C. **Do not freeze.**
2. RPR antigenic suspension contains thimerosal at 0.095% as preservative, which is toxic if ingested.
3. The source material used in the preparation of the controls has been tested and confirmed negative for HCV, HIV and HBsAg. Because no test can offer complete assurance that products derived from human source will not transmit infectious agents, it is recommended that these controls be treated as potential biohazard materials.
4. 0.095% Sodium azide may react with lead and copper plumbing to form highly explosive salts. On disposal, flush with large quantities of water.
5. Disposal in accordance with local regulations.
6. Follow the good laboratory practices (GLP) related to storage, handling and material disposal.

WARRANT

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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