MS 10310030096

Imuno-Rápido Dengue IgG/IgM

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ENGLISH

SUMMARY

Dengue viruses is transmitted by the mosquito Aedes aegypti. There are four known distinct serotypes (Den1, Den-2, Den-3 and Den-4). Incubation period is

from 3 to 15 days, the average is from 5 to 6 days. Dengue infection is associated with high fever (39° to 40°C), headache, muscle pain, myalgia, morbidity, arthralgia, retro-orbital pain, asthenia, anorexia, nausea, puke, exanthema and skin rash. Dengue hemorrhagic fever is the most serious disease and includes gingival hemorrhage, petechia and ecchymosis, gastroenterorrhagic, shock and death. The most Laboratory finding is thrombocytopenia. In 2 to 4% of the patients infected a second time, dengue hemorrhagic fever is more likely to occur.

Dengue is considered to be the most important viral disease due to the human morbidity and mortality it causes. Dengue is widely distributed throughout the tropical and subtropical areas of the world. Epidemic period is likely to occur during or immediately after the rainy seasons. It is known that around 2.5 millions of people live in areas where the transmission of Dengue is more likely to occur. In Brazil, from 1990 to 2000, Den-1 and Den-2 serotypes disseminated mainly in urban areas of southeast and northeast. In 2003, Den-1, Den-2 and Den-3 serotypes were researched in 23 Brazilian states and was registered a crescent number of Dengue hemorrhagic fever. Dengue is endemic, with epidemic peaks in the first four months of the year. The inefficient vector control associated with the lack of an efficient vaccine increased the causes of Dengue hemorrhagic fever.

Nowadays, rapid test is the safest form to diagnose the primary and secondary infection of Dengue and it is essential for patient management. Such results are associated to immune response, which includes IgM antibodies produced by 5° to 8° day of the symptoms and it persists for 30-60 days, and indicates a recent infection or acute phase of the disease.

IgG appears on the 14° day and it can persist for life. As the infection for one subtype does not give immunity to other subtypes, secondary infections show that IgG rises within 1 to 2 days after the onset of symptoms and induces IgM specifier later.

Wama Diagnostica Imuno-Rápido DENGUE IgG/IgM is a solid phase immunochromatographic assay for the rapid, qualitative and differential detection of IgG and IgM antibodies against the four serotypes of Dengue virus in human serum and plasma (EDTA, Sodium citrate and Heparin).

METHOD PRINCIPLE

Specific immunoglobulin of Dengue IgG and/or IgM present in the sample bind to recombinant antigens (Den-1, Den-2, Den-3 and Den4) of viral envelope conjugated to colloidal gold. This reaction forms a complex antigen-antibody which migrates along the length of the test device by capillary action and is captured by anti-human IgG and/or anti-human IgM immobilized in two test lines and generates a pink colored band.

If IgG and/or IgM antibodies to Dengue virus are not present in the sample, there is no color appearance in the respective areas. A reagent control immobilized on the membrane of the test device will determine the appearance of a third pink colored band in case the test procedure is performed properly and the test reagents of control line are working.

KIT PRESENTATION

REF 626010-R: (10 determinations)

1. Test Device : 10 units
2. Assay Diluent 1x 1.2ml
3. Capillary Pipette : 10 units
4. Package Insert

REF 626020-R: (20 determinations)

1. Test Device : 20 units
2. Assay Diluent 1x 3.0ml
3. Capillary Pipette : 20 units
4. Package Insert

REF 626025-R: (25 determinations)

1. Test Device : 25 units
2. Assay Diluent 1x 3.0ml
3. Capillary Pipette : 25 units
4. Package Insert

REF 626040-R: (40 determinations)

1. Test Device : 40 units
2. Assay Diluent 2x 3.0ml
3. Capillary Pipette : 40 units
4. Package Insert

REAGENT STABILITY

Test Device (1): Ready for use. Stable if stored at 2-30°C up to expiration date. **Assay Diluent (2):** Phosphate buffer. Contains Sodium azide 0.095%. Ready for use. Stable if stored at 2-30°C up to expiration date.

Capillary pipette (3): 5µl Plastic capillary pipette. **Obs.:** The kit presents good performance after being used for the first time. It's stable up to the expiration date if stored at 2-8°C.

MATERIAL REQUIRED BUT NOT PROVIDED

- Automatic pipette
- Paper towel
- Discard material container
- Paper towels

SPECIMEN COLLECTION AND HANDLING

Serum or plasma specimens may be used with this test. Plasma obtained from EDTA, Sodium citrate and Heparin and free of haemolysis, lipaemic and contamination are preferable. Fresh specimens are recommended, otherwise, store the specimens refrigerated at 2-8°C for 1 week. If longer storage is required, the specimens should be stored in freezer at -20°C.

ATTENTION

If the specimen was frozen, it should be completely thawed, homogenized and placed in a vertical position to allow the sedimentation of any particle. Do not shake the specimen. Diluted specimen can cause false negative result.

PROCEDURE

- Let test device, specimens and assay diluent to reach to room temperature, if they are stored refrigerated, prior to test.
- Remove test device from individually sealed pouch and use it immediately.
- Place test device on a dry and flat surface.
- Using a micropipette or a capillary pipette (3), add 5µl of serum or plasma into the sample well marked "S".

Obs. When using capillary pipette (3), squeeze the pipette and place it into the specimen, let the specimen reach the mark on the pipette slowly. Slightly touch the pipette on sample well marked "S" and squeeze it to dispense the specimen on membrane.

- Add 3 drops (90µl) of assay diluent into the assay diluent well round shaped.
- Interpret test results from 15 to 20 minutes. Do not read the results after 20 minutes.

INTERPRETATION

NEGATIVE

One pink band appears in the control area (C). No Dengue infection was detected. Make a new test after 3 to 5 days if dengue infection is apparent.

IgM POSITIVE (Primary dengue infection)

Two pink bands appear, one in the (M) area and the other in the control area (C).

IgG POSITIVE (Secondary or past dengue infection)

Two pink bands appear, one in the (G) area and the other in the control area(C).

IgG and IgM POSITIVE (Late primary or early secondary dengue infection)

Three pink bands appear in the (M), (G) and (C) areas.

OBS. Any colored band in the (M) and (G) areas should be interpreted as POSITIVE.

INVALID No colored band appears in the (M), (G) and (C) or no band in the control area (C). Do not take into account any result after the reading time recommended.

EXPECTED VALUE

Imuno-Rápido DENGUE IgG/IgM was compared to inhibition test of hemagglutination and the results are showed below:

SENSITIVITY AND SPECIFICITY

	102 IgG/IgM positive specimens		
	Positive	Negative	Sensitivity
IMUNO-Rápido Dengue IgG/IgM	101	1	101/102x100%=99%

	200 IgG/IgM negative specimens		
	Positive	Negative	Sensitivity
IMUNO-Rápido Dengue IgG/IgM	4	196	196/200x100%=98%

Sensitivity: 99%

Specificity: 98%

PRECISION

Intra-assay

Within-run precisions have been determined by the testing of four specimens 10 times: 1 negative, 1 low positive, 1 medium positive and 1 high positive. All values were correctly identified 99% of the time.

Inter-assay

Between-run precisions were determined by 10 independent testing on the same 4 specimens: 1 negative, 1 low positive, 1 medium positive and 1 high positive.

Three different lots of IMUNO-Rápido DENGUE IgG/IgM were tested for 3 months using negative, low positive, medium positive and high positive specimens. All values were correctly identified 99% of the time.

LIMITATION OF THE TEST

1. IMUNO-Rápido DENGUE IgG/IgM (human serum and plasma) is for in vitro diagnostic use only. This test detects the presence of antibodies anti-Dengue and it should be used only with human serum and plasma. Neither quantitative value nor antibodies anti-Dengue concentrations can be determined since this test is qualitative.

2. IMUNO-Rápido DENGUE IgG/IgM detects the presence of antibodies anti-Dengue in the specimens and should not be used as the sole criteria for the diagnosis of Dengue virus infection.

3. As with all diagnostic tests, all results must be considered with other clinical information available to the physician.

4. If the test is negative and clinical symptoms persist, additional follow-up using other clinical methods is recommended. A negative result does not preclude the possibility of infection of dengue virus.

PRECAUTIONS AND WARNINGS

- For in vitro diagnostic use only.
- The instructions must be read prior to testing.
- The expiration date is indicated by the last day of the month printed on the test device pouch and on the kit box. Do not use kit components beyond the expiration date.
- Protect the kit from heat and sunlight.
- Do not freeze the test device.
- Sodium azide 0.095% may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquid, flush with large volumes of water to prevent azide buildup. Sodium azide might be toxic if ingested.
- Allow all kit components and specimen reach to room temperature prior to testing.

8. Do not interchange kit components with those from sources other than the same catalog number from WAMA DIAGNOSTICA.

9. Disposal in accordance with local regulation.

10. Follow good laboratory practices (GLP) in storing, dispensing and disposing of these materials.

WARRANT

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama´s technical support. The warranty will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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