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ENGLISH

CLINICAL IMPORTANCE

Prostate Specific Antigen (PSA) is produced by the prostate gland and by epithelial cells. It is an intracellular glycoprotein (M.W. 34,000 dalton). When abnormal events occur in the prostate, it may result in high levels of PSA in the blood which help the diagnosis of prostate cancer and benign prostatic hyperplasia. Inflammatory process such as prostatitis also changes PSA levels which became normal after the treatment.

It is established that the PSA levels in the blood rise according to the age of patient, thus, dependent age. A major problem on PSA result interpretation is related to the CUT-off value to be used in the prostate cancer diagnosis. There are several recommendations for a better result interpretation such as rectal exam and the rise of blood levels of PSA in a short period of time.

Many authors have demonstrated prostate cancer detection by biopsy in 25% of cases where present PSA levels from 2,5 to 4,0 ng/ml, 30% from 4 to 10ng/ml and 65% above 10ng/ml, men above 50 years of age. As blood levels of PSA rise with the age, cut-off of 2,5ng/ml should be always regarded as individuals below 50 years of age.

Imuno-Rápido PSA was developed to assist the actual recommendations, with sensitivity of 2,5ng/ml.

It is recommended that positive results be confirmed against other methodologies.

PRINCIPLE OF THE METHOD

PSA existent in the specimen binds to anti-PSA antibody conjugated to colloidal gold. This mixture migrates by capillary action through the membrane and reacts with anti-PSA antibodies fixated in the membrane and generates a colored line in the test region "T".

A control reagent immobilized in the membrane will determine the appearance of a second colored line in the control region "C". This reaction will demonstrate that reagents are working properly.

KIT PRESENTATION

REF 630010-R - 10 determinations

1. Test device: 10 units
2. Diluent solution (1x0.5ml)
3. Instructions for use

REF 630020-R - 20 determinations

1. Test device: 20 units
2. Diluent solution (2x0.5ml)
3. Instructions for use

REF 630040-R - 40 determinations

1. Test device: 40 units
2. Diluent solution (4 x 0.5ml)
3. Instructions for use

PREPARATION AND STABILITY OF REAGENTS

TEST DEVICE (1): Ready for use. Stable up to expiration date if stored at 2–30°C. **Do not freeze.**

DILUENT SOLUTION (2): Ready for use. It contains sodium azide 0.095%. Stable up to expiration date if stored at 2–30°C. **Do not freeze.** The solution should be brought to room temperature prior to use.

Obs.: Kit presents consistent results after its first handling and it is stable up to expiration date if stored at 2-30°C.

SPECIMEN COLLECTION AND STORAGE

IMUNO-RÁPIDO PSA may be performed using serum.

The test should be performed immediately after specimen collection. It is not recommended to keep the specimens at room temperature for prolonged periods. The serum specimens may be stored at 2 to 8°C for up to 3 days. For prolonged storage, the specimens should be stored at 2-20°C.

It is recommended to avoid the use of icteric specimens since they can present false positive results.

The specimens should be brought to room temperature prior to use. Freeze specimens should be completely thawed and homogenized prior to use.

Avoid repeatedly thawing and freezing the specimens.

MATERIAL NEEDED BUT NOT PROVIDED

- Container for specimen collection
- Automatic pipette
- Paper towel
- Container for material disposal

PROCEDURE

All components should be brought to room temperature prior to use.

1. Test device should be brought to room temperature prior to remove it from the foil pouch.
2. Pipette 50µl of specimen (no air bubbles) into the specimen well (>) of the test device
3. Add 1 drop of diluent solution (50µl)
4. Interpret the test results at 10 minutes. Do not interpret test results


after 15 minutes.

Obs.: Specimens ranging from 1.8 to 2.2ng/ml can present results weakly positive (shadow), with interpretation after 10 minutes.

INTERPRETATION OF THE TEST


NEGATIVE:

The presence of only one colored band in the control region (C).


POSITIVE:

The presence of two colored bands, one in the control region (C) and the other in the test region (T).

INVALID:

If no colored band appears in the test region (T) and control region (C), or if no colored band appears in the control region (C). It is recommended that the specimen be re-tested.

Insufficient specimen volume or the instructions may not have been followed correctly. It is recommended to verify the procedures and to re-test using a new test device.

Do not interpret the test results after the prescribed period.

QUALITY CONTROL

A control procedure is included in the test. A colored band which appears in the control region (C) is regarded as an internal control procedure. It confirms that the specimen volume is sufficient and the instructions are followed correctly.

The standard controls are not provided. It is recommended to test positive and negative controls as good manufacturing practices in order to confirm the procedure and test performance.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity:

Analytical sensitivity of 2.5ng/ml. (lower value detected).

Clinical or diagnostic sensitivity:

250 specimens were tested, 50 positive specimens (above 2.5ng/ml). The kit presented a false negative result which determines a sensitivity of 98%.

Specificity:

250 specimens were tested, 200 negative specimens (0.005 – 2.4ng/ml). All specimens were correctly identified with no false positive which determines a specificity of 100%.

Precision

Within run

The repeatability was determined by using pilot lots in three assays under the same conditions and using 10 specimens of known results in triplicate. The specimens were correctly identified 100% of the time.

Between-run

The reproducibility was determined by performing tests with pilot tests in three different days with 10 specimens. The specimens were identified 100% of the time.

LIMITATIONS OF THE TEST

IMUNO-RÁPIDO PSA is for *in vitro* diagnostic use only. The test may be performed using to detect PSA in serum.

PSA test should not be used as a unique criterion for the diagnosis of prostate cancer.

A significant number of patients with BPH, more than 15% and lower than 1% of healthy individuals has elevated PSA concentrations. Other clinically available tests are required if questionable results are obtained and should only be made by the physician after all clinical and laboratory findings have been evaluated.

PSA concentrations are not trustful in patients submitted to hormone therapy or prostate gland manipulation.

PRECAUTIONS AND WARNINGS

1. For "*in vitro*" diagnostic use only.
2. Read carefully the instructions for use prior to use.
3. Expiration date is the last day of the month printed on labels.
4. The test kit should be kept away from direct sunlight
5. Do not freeze test device.
6. The kit contains 0.095% sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive salts. On disposal, flush with large quantities of water. Sodium azide may be toxic if ingested.
7. Do not use kit beyond the expiration date.
8. Do not interchange kit components from different kits.

9. Disposal in accordance with local regulations.

10. Follow good laboratory practices (GLP) related to storage, handling and material disposal.

WARRANT

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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1. Catalona WJ, Southurick PC, Slawin KM, et al., Comparison of percent free PSA, PSA density and age-specific PSA cut-offs for prostate cancer detection and staging. Urology 1:56(2):255-60, 2000.
2. Christens A, Laurell CB, Lija H. Enzymatic activity of prostate – specific antigen and its reaction with extracellular serine proteinase inhibitors. Eur J Biochem; 194:755-763, 2000
3. Lodding P, Aus G, Bergdahl S, et al. Characteristics of screening detected prostate cancer in men 50 of 66 years old with 3 to 4ng/mL prostate specific antigen. Journal of Urology 159(3):899-903, 1998
4. Wang MC, Valenzuela LA, Murphy GP, et al., Purification of human prostate specific antigen. Invest Urol; 17:159-163, 1979.

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