MS 10310030084

Imuno-Rápido

SANGUE OCULTO

Dispensa Dieta Alimentar

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## ENGLISH

### SUMMARY

The presence of hemoglobin in feces is indicative of early gastrointestinal cancer and it is frequently found in patients with anaemia caused by the lack of iron.

Therefore, the detection of blood in feces is an important phase of diagnosis. Abnormal bleedings can be present from stomach to rectum.

Colorectal cancer is the third most common cancer in the USA, approximately 135,000 cases a year, and it is the second in cancer death. Many cancers, if not all of them, develop from a benign adenomatous polyps for more than 10 years, although approximately 5% of the polyps become malign.

Fecal Occult Blood test is a valuable method to detect gastrointestinal cancer and adenomatous polyps through the identification of patients for further diagnosis (colonoscopy).

There are many methods to detect occult blood in feces.

**Imuno-RÁPIDO SANGUE OCULTO** is an immunochromatographic test which employs a combination of a marked monoclonal antibody and a polyclonal antibody anti-hemoglobin human of solid phase using a high degree of sensitivity and specificity. No food restrictions to specimen collection.

### PRINCIPLE OF THE METHOD

Hemoglobin present in the specimen binds to the conjugate monoclonal antibody forming an antigen-antibody complex. The complex migrates on the membrane chromatographically by capillary action to react with anti-hemoglobin antibody human in the test area (T) and produces a colored line. The absence of hemoglobin, will mean no colored line in the test area (T).

The mixture migrates to the control area (C) of the membrane. The conjugate which is not binded to antigen joins to the reagents of this area producing a light pink colored line. This process confirms that the reagents and device are functioning correctly.

### KIT PRESENTATION

REF 623010-R (10 determinations)

1. Test device: 10 units

2. Specimen collection tube: 10 units

3. Instructions for use

REF 623020-R (20 determinations)

1. Test device: 20 units

2. Specimen collection tube: 20 units

3. Instructions for use

REF 623040-R (40 determinations)

1. Test device: 40 units

2. Specimen collection tube: 40 units

3. Instructions for use

REF 623080-R (80 determinations)

1. Test device: 80 units

2. Specimen collection tube: 80 units

3. Instructions for use

REF 623010T-R (10 determinations)

1. Strip-Test: 10 Units

2. Collectors Sample (2ml): 10 Units

3. Card-Test: 10 Units

4. Instruction for Use

REF 623020T-R (20 determinations)

1. Strip-Test: 20 Units

2. Collectors Sample (2ml): 20 Units

3. Card-Test: 20 Units

4. Instruction for Use

REF 623040T-R (40 determinations)

1. Strip-Test: 40 Units

2. Collectors Sample (2ml): 40 Units

3. Card-Test: 40 Units

4. Instruction for Use

### MATERIAL REQUIRED BUT NOT PROVIDED

• Specimen collection container

• Disposable material container

### PREPARATION AND STABILITY REAGENT

• **TEST DEVICE (1):** Keep at 2-30°C. **Do not freeze.**

• **SPECIMEN COLLECTION (2):** It contains extraction buffer. Keep at 2-30°C. **Do not freeze.**

The kit presents good performance after being used for the first time. It is stable up to the expiration date if store at 2-30°C.

### SPECIMEN COLLECTION

Collect the specimen of feces from one or three defecations (preferable specimen).

**Specimen collection procedure: (Collected by the patient at home or by a specialist at laboratory):**

1. Poke the stick of the Collection Tube (2) into the fecal specimen three times at different sites.

2. Reinsert the stick in the tube and tighten securely.

3. Store the Specimen collection at 4-8°C (when collected at home)

**OBSERVATION:** If the specimen is collected in three days consecutively (one defecation a day), the procedure above should be repeated with each specimen using the same collector.

4. Send the Specimen Collection (2) to laboratory up to 8 days (when collected at home).

### PROCEDURE

1. Remove the Test Device (1) from the sealed foil pouch.

2. Break the tip of the Specimen Collection and transfer 2 drops of the specimen in the (S) area of the test device.
3. Read result between 10 to 15 minutes.


#### Test Strip

1. Allow the test strip reach room temperature before removing it from the laminated pouche.
2. Peel off adhesive and stick the test strip in the middle of the car with the arrows pointing downward.
3. Break the tip of Specimen collection and transfer 2 drops of the specimen on the base of the test strip
4. read result between 10 to 15 minutes.

#### READING

#### NEGATIVE:


One colored line is visible in the control area ©.

#### POSITIVE:

Two colored lines are visible, one in test area (T) and another in control area ©.

Any shade of pink line in the test area indicates positive result.

**INVALID:** If no line is visible in test area (T) and in control area ©, or if no line is visible in control area ©. The test should be repeated.

#### TEST PERFORMANCE

Analytic sensitivity=50ng/ml of hemoglobin

It was evaluated 135 symptomatic patients undergone colonoscopy and presenting some types of alteration, with possible lack of blood in the feces (polyps >0.5cm, colorectal cancer, haemorrhoids, Crohn's disease, anal fissure, ulcer)

**Imuno-RÁPIDO SANGUE OCULTO** and HEMOCULT II (SKD, California, USA) were evaluated and compared and the sensitivity and specificity are as follow:

**Specificity:** When the criteria was bleeding for several pathologies, **Imuno-RÁPIDO SANGUE OCULTO** specificity was 94% against 88% of HEMOCULT II. When the criteria was only the colorectal cancer and polyps >0.5cm, the specificity was 91% and 80% respectively

**Sensitivity:** When the criteria was bleeding for several pathologies, **Imuno-RÁPIDO SANGUE OCULTO** sensitivity was 38% against 26% of HEMOCULT II.

When the criteria was only the colorectal cancer and polyps >0.5cm, the sensitivity was 52% and 28% respectively.

**Cross reaction:** It was not observed cross reactions with hemoglobin of bovine, rabbit, horse and sheep.

**Work Zone:** 0.04mg/gr to 120mg/gr of feces was determined as work zone.

Specimens collected with formalin buffered preservative are not appropriated for recal occult blood rapid test. This test should only be performance using fresh specimen of feces.

#### LIMITATION OF USE

The presence of blood in the feces can present some etiologies besides colorectal bleeding such as haemorrhoids, gastric irritation etc.

**Imuno-Rápido SANGUE OCULTO** detects human blood in quantities bigger than physiological loss of the healthy individuals.

Negative results does not exclude bleeding since it can be intermittent.

Early phase of colorectal polyposis can show no bleeding. For this reason, it is recommended to annually perform the test in individuals above 45 years.

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be interpreted with other pieces of information available to the physician and after clinical, radiologic and laboratorial findings have been evaluated.

Specimens collected with formalin buffered preservative are not appropriated for fecal occult blood rapid test. This test should only be performed using fresh specimens of feces.

#### PRECAUTIONS AND WARNINGS

1. For in vitro diagnostic use only.
2. Follow all steps on the package insert prior to use the kit.
3. Specimen collection can irritate the skin, eyes and mucous. In case of contact, wash with abundant water
4. The expiration date matches the last day of the month printed on the label on test device and kit box.
5. Avoid keeping the kit at high temperatures and direct sunlight.
6. Do not freeze the test device.
7. Do not use the kit beyond the expiration date.
8. Disposal in accordance with local regulations.

9. Follow Good Manufacturing Practices (GMP) at storing, handling and material disposable.

#### WARRANT

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replace if technical support finds evidence that running, handling and storage were not properly followed.

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#### BIBLIOGRAFIA / BIBLIOGRAPHY / BIBLIOGRAFIA

1. Hardcastle, J. D. et al.: Randomized controlled trial of faecal-Occult screening for colorectal cancer. **Lancet**, 348: 1472-1477, 1996.
2. Helm, J. F. and Sandler, R. S. Screening do cáncer colorretal. In: Lang, R. S. and Isaacson, D. H. (Editors). **Clínicas Médicas da América do Norte**, 83 (6): 1279-1296, 1999.
3. Kronborg, O. et al.: Randomized study of screening for colorectal cancer with faecal-occult-blood test. **Lancet**, 348: 1467-1471, 1996.
4. Mandel, J. S. et al.: Sensitivity, specificity, and positive predictivity of the Hemmocult test in screening for colorectal cancers. **Gastroenterology**, 97:597-600, 1989.
5. Mandel, J. S. et al.: Reducing mortality from colorectal cancer by screening for faecal occult. **N. Engl. J. Med.**, 328: 1365-1371, 1993.
6. Price, C. P. et al.: Disposable integrated immunoassay devices. In: Price, C. P. and Newman, D. J. (Editors). **Principles and Practice of Immunoassay, 2nd ed., Macmillan Reference**: 581-603, 1997.
7. Yong, G. P. et al.: Faecal occult blood tests: choice, usage and clinical applications. **Clin. Biochem. Revs.** 13:161-169, 1992.
8. Winawer, S. J. et al.: Screening for colorectal cancer with fecal occult blood testing and sigmoidoscopy. **J. Natl. Cancer Inst.** 85: 1311-1318, 1993.

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