MS 10310030083

# Imuno-Rápido CHLAMYDIA

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## ENGLISH

### SUMMARY

Chlamydiae are among the most common pathogens of animal kingdom. They are gram-negative bacteria which are adapted to obligate intracellular existence in eukaryotic host cells. Chlamydiae replicate in cytoplasm of host cells by forming intracellular inclusions and differ from the virus since they have DNA and RNA.

They belong to Chlamydiaceae family, with a genus, *Chlamydia*, which contains four species: *C. trachomatis*, *C. psittaci*, *C. pecorum* and *C. pneumoniae*.

*C. trachomatis* is naturally found living only inside human cells and it is one of the most common sexually transmitted infections in people worldwide.

It is responsible for 50% of non-gonococcal urethritis and an even larger percentage among post-gonococcal urethritis cases. Chlamydia can also cause proctitis in homosexual, conjunctivitis (associated with genital infection through genital-manual-ocular contact), venereal lymphogranuloma and asymptomatic infection. Cervicitis, urethritis, endometritis, salpingitis,

proctitis, conjunctivitis, venereal lymphogranuloma, and asymptomatic infection are more common in women.

Salpingitis can cause infertility and ectopic pregnancy. Endocervical infection contaminates 60% of newborns through vaginal parturition causing conjunctivitis (18% to 50%) and pneumonia (3% to 15%) on infants.

There are several available methods for *Chlamydia* infections diagnosis such as the test using colored specimen, cell culture, detection of antigens using antibodies labelled with FITC (immunofluorescence) or enzyme (immunoenzymatic), molecular (PCR) and antibody detection.

Imuno-RÁPIDO CHLAMYDIA is an immunochromatographic test, antigen detection, using specific antibodies which selectively identify *Chlamydia trachomatis* with high degree of sensitivity and specificity.

### PRINCIPLE OF THE METHOD

After the specimen collection, if *Chlamydia trachomatis* is present, the antigen binds to a dye antibody conjugate (gold colloid) by forming an antigen-antibody complex. This complex migrates through test card and binds to an antibody in the test region (T), a light pink band appears. On the absence of *Chlamydia trachomatis*, no colored band will appear in the test region (T).

The mixture migrates to the control region (C). The conjugate, not bound to the antigen, is mixed with reagents of (C) region and a light pink band appears showing that the reagents are properly running.

### KIT PRESENTATION

**REF 618010-R (10 determinations)**

1. Test device: 10 units
2. Swabs: 10 units
3. Reagent A: 1x4ml
4. Reagent B: 1x4ml
5. Pipette
6. Test tubes with dropper tips: 10 units
7. Workstation
8. Instructions for use.

**REF 618020-R (20 determinations)**

1. Test device: 20 units
2. Swabs: 20 units
3. Reagent A: 2x4ml
4. Reagent B: 2x4ml
5. Pipette
6. Test tube with dropper tips: 20 units
7. Workstation
8. Instructions for use.

**REF 618040-R (40 determinations)**

1. Test device: 40 units
2. Swabs: 40 units
3. Reagent A: 3x4ml
4. Reagent B: 3x4ml
5. Pipette
6. Test tube with dropper tips: 40 units
7. Workstation
8. Instructions for use.

### MATERIAL REQUIRED BUT NOT PROVIDED

- Swab tubes.
- Swab or cytologic brush (urethral specimens)
- Centrifuge.
- Disposable material container.

### REAGENT STABILITY AND STORAGE

• TEST DEVICE, REAGENT A e B: They should be kept at 2-30°C. Stable up to expiration date. Do not freeze.

### SPECIMEN COLLECTION AND STORAGE

The quality of specimen obtained is extremely important for the accuracy of the results of any diagnostic test for *Chlamydia trachomatis*.

*Chlamydia trachomatis* is an intracellular parasite which grows inside the epithelial cells of cervix and urethra. For an adequate specimen, firmly rotate the swab against the urethral or cervical canal. Therefore, the specimen collection methods below are recommended:

#### CERVICAL SPECIMENS:

- 1) Use the swab provided on the kit. Alternatively, any fiber-tip swab may be used. Do not use wood-shaft swab, cotton-tip swab or swabs that contain agar or gelatin as means of transport.
- 2) Remove the excess of mucus from the endocervical area with a cotton-tip swab and discard it. The excess of mucus can cause false positive results.
- 3) The swab should be inserted into the endocervical canal. Firmly rotate the swab against the canal for 15 to 30 seconds.
- 4) Then, remove the swab with no vaginal contact. If the test is to be conducted immediately, put the swab into the extraction tube.

#### MALE URETHRAL SWAB SPECIMENS:

Note: Instruct patients not to urinate for at least one hour prior to specimen collection.  
1) Fiber-tip swab or cytologic brush (not provided) should be used. Do not use swab with agar and gelatin as means of transport.

- 2) Insert the swab into the urethra about 2-4cm, rotate it for 10 seconds.
- 3) Then, remove the swab. If the test is to be conducted immediately, put the swab into the extraction tube.

#### MALE URINE SPECIMENS:

1) Instruct patients to collect 15-30ml of clean first morning urine in a sterile urine cup with no preservative. First morning urine specimens are preferred to achieve the highest concentrations of *Chlamydia* antigen, in case of positive specimens.

Observations: The excess of salt (uric acid) in urine, evidenced by a large amount of precipitates, can interfere on test results.

#### Observations:

It is recommended that all specimens are tested as soon as possible after collected. If immediate testing is not possible, the swab should be placed in a dry transport tube for storage. The swab may be stored for 4-6 hours at room temperature (15-30°C) or 24-72 hours refrigerated (2-8°C). The urine specimens can be stored refrigerated (2-8°C) for 24 hours. Do not freeze. All specimens should be allowed to reach room temperature before testing.

#### PROCEDURE

##### I. SPECIMEN EXTRACTION TO SWAB:

1. Add 5 drops of reagent A (approximately 300µl) to the extraction tube. This reagent is colorless.
2. Insert the swab into the extraction tube and mix vigorously for 2 minutes. Let stand for 2 minutes.
3. Fill in the pipette for reagent B up to the marked line (approximately 200µl) then add the

reagent B to the extraction tube, this reagent is cloudy. Mix vigorously the swab for 1 minute, the solution will present a blue or greenish color. If the swab is bloody, the color will turn to yellow or orange. Blood excess can cause false positive result. Let stand for 1 minute.

4. Press the swab against the side of the tube to remove the excess of solution and remove the swab.

5. Fit the dropper tip on top of the extraction tube.

### III. SPECIMEN EXTRACTION FOR MALE URINE:

1. Mix gently the urine. Transfer 5ml of urine to a test tube and add 5ml of distilled water. Centrifuge for 15 to 20 minutes at 3,000 rpm.

2. Carefully discard the supernatant and dry the rim of tube by blotting onto absorbent paper. If the supernatant is not discarded, it can cause false positive results.

3. Fill in the pipette for reagent B up to the marked line (approximately 200µl) then add the reagent B to test tube with urine pellet. Mix vigorously with the pipette up to a complete mixture.

4. Transfer all the mixture to the extraction tube and let stand for 1 minute

5. Add 5 drops of Reagent A (approximately 300µl) to the extraction tube. Mix vigorously for 2 minutes and let stand for 2 minutes more.

6. Fit the dropper tip on top of the extraction tube.

### III. TEST


1. Bring the test device to room temperature prior to testing.

2. Add 3 drops (approximately 100µl) of extracted solution in the specimen region (S) of test device.

3. Interpret the results after 10 minutes. Do not interpret them after 20 minutes.


### INTERPRETATION OF THE RESULTS

**NEGATIVE:**  
One colored line appears in the Control region (C).



### POSITIVE:

Two colored lines appear. One line should be in the control region (C) and another line should be in the test region (T). The intensity of color in the test region (T) may vary, but it should be considered positive whenever there is even a faint colored line

### INVALID:

No band appears in the test region (T) and in the control region (C), or no band appears in the control region (C). Review the procedure and repeat the test with a new specimen or with the same extracted solution.

### TEST PERFORMANCE

Immuno-RÁPIDO CHLAMYDIA is a screening qualitative test for the detection of *Chlamydia trachomatis*.

### SENSITIVITY AND SPECIFICITY:

Immuno-Rápido Chlamydia was compared to PCR using negative and positive specimens. The results show that the Immuno-Rápido Chlamydia from Wama Diagnostica presented high concordance with PCR assay, indicating that this test has a higher specificity and sensitivity, as the results showed below:

Endocervical specimens:  
Sensitivity: 88,5% - Specificity: 96,7%

Male urethral specimens:  
Sensitivity: 78,4% - Specificity: 92,9%

Male urine specimens:  
Sensitivity: 90,9% - Specificity: 99,0%

### CROSS REACTIVITY:

The following organisms, in concentration of 10<sup>6</sup> cells/ml, were found negative when tested with Immuno-Rápido Chlamydia.

- *Neisseria gonorrhoea*
- *Neisseria meningitidis*
- *Proteus vulgaris*
- *Candida albicans*
- *Pseudomonas aeruginosa*
- *Acinetobacter calcoaceticus*
- *Staphylococcus aureus*
- *Enterococcus faecium*
- *Gardnerella vaginalis*
- *Branhamella catarrhalis*
- *Enterobacter aerogenes*
- *Proteus mirabilis*
- *Streptococcus faecalis*
- *Klebsiella pneumoniae*
- *Salmonella choleraesuis*
- *Hemophilus influenzae*

### LIMITATIONS OF USE

Immuno-RÁPIDO CHLAMYDIA should be used for the detection of *Chlamydia trachomatis* antigen from endocervical and urethral areas and from male urine specimens. Other specimens collected from rectum and conjunctive are not recommended.

As with all immunological tests for *Chlamydia trachomatis*, this test does not differentiate active and inactive organisms.

Excessive mucus or blood on the swab may cause false positive results, therefore, specimen collection is extremely important for the test performance.

The test results should be interpreted in conjunction with other laboratory and clinical data available to the physician.

Other tests should be done such as cell culture when a negative result is found with Immuno-Rápido Chlamydia on the presence of persistent clinical symptoms.

### PRECAUTIONS AND WARNINGS

1. In vitro diagnostic use only.

2. Read carefully the package insert prior to using the test.
- 3. Reagents A and B may cause skin, eyes and mucus irritation. Wash the parts in contact with the reagents with abundant water.
- 4. Expiration date corresponds to the last day of the last month printed on the pouche of the test device and on the kit box.
- 5. Keep the kit away from the sunlight and heat.
- 6. Do not freeze the test device.
- 7. Bring the reagents to room temperature (20-25°C) prior to using the test.
- 8. Do not use the kit after the expiration date.
- 9. Disposal in accordance with local regulations.
- 10. Follow the good manufacturing practices (GMP) on storing, handling and material discard.

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### BIBLIOGRAPHY

1. Barnes, R. C.: Laboratory diagnosis of human chlamydial infections. Clin. Microbiol. Rev., 2: 119, 1989.
2. Brunham, R. C. et al.: *Chlamydia trachomatis* infection in women with ectopic pregnancy. Obstet. Gynecol., 67: 722, 1986.
3. Jaschek, G. et al. Direct detection of *Chlamydia trachomatis* in urine specimens from symptomatic and asymptomatic men by using a rapid polymerase chain reaction assay. J. Clinical Microbiology, 31: 1209-1212, 1993.
4. Martin, D. H.: Infeções por *Chlamydia*. Clin. Med. Am. Norte, 6: 1479-1500, 1990.
5. Pearlman, M. D. and McNeely, S. G.: A review of the microbiology, immunology, and clinical implications of *Chlamydia trachomatis* infections. Obstet Gynecol. Surv., 47: 448, 1992.
6. Price, C. P. et al.: Disposable integrated immunoassay devices. In: Price, C. P. and Newman, D. J. (Editors). Principles and Practice of Immunoassay, 2nd ed., Macmillan Reference: 581-603, 1997.
7. Sanders JW, et al. Evaluation of an enzyme immunoassay for detection of *Chlamydia trachomatis* in urine of asymptomatic men. J. Clinical Microbiology, 32: 24-27, 1994.
8. Schachter J. Sexually transmitted *Chlamydia trachomatis* infection. Management of the most common venereal disease. Postgraduate Medicine, 72: 60-69, 1982.
9. Schachter, J. and Stamm, W.E.: *Chlamydia*. In: Murray, P.R.; Baron, E. J.; Pfaller, M. A.; Tenover, F. C. and Tenover, R. H. (editors). Manual of Clinical Microbiology, 7th ed. ASM Press, Washington, D. C.: 795-806, 1999.
10. Zhang, J. P. and Stephens, R. S.: Mechanism of *C. trachomatis* attachment to eukaryotic host cells. Cell, 69: 861-869, 1992.

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