MS 10310030091

Uri-Color Check

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ENGLISH

SUMMARY

Urine test assay aids the physician on accurate pieces of information on renal pathology and some extra-renal diseases as well. This test has been used for several years due to its easy-to-use procedure, low cost and abundant sensitivity.

URI-COLOR CHECK from WAMA Diagnóstica detects accurately semi-quantitatively the presence of: urobilígen, glucose, bilirubin, ketones, pH, blood/hemoglobin, Specific Gravity, protein, nitrite and leucocytes in urine specimen.

PRINCIPLE OF THE METHOD

Urobilígen: The test is based on the diazotization reaction of 4-methoxybenzene diazonium salt and urinary urobilígen in a strong acid medium. The color changes range from pink to brown-red. **Glucose:** This test is based on a sequential enzyme reaction. First, glucose oxidase catalyses the formation of gluconic acid and hydrogen peroxide from the oxidation of glucose. A second enzyme, peroxidase, catalyses the reaction of hydrogen peroxide with potassium iodide chromogen. Colors change range from blue through greenish-brown and brown to dark-brown.

Ketones: This test is based on the reaction of acetoacetic acid in the urine with nitroprusside. The resulting color ranges from tan when no reaction takes place and to purple for positive reaction. **Bilirubin:** This test is based on the coupling of bilirubin with 2,4-dichlorobenzene diazonium Na in a strong acid medium. The color changes from light tan to brown.

Protein: This test is based on the color changing of the indicator, tetramotofenol, in the presence of protein. A positive reaction is indicated by a change on the color from yellow to green and then to greenish-blue.

Nitrite: This test is based on the reaction of p-arsanilic acid and nitrite (which is derived from dietary nitrate in the presence of bacteria) in urine to form a diazonium compound. The diazonium compound, in turn, couples with N-(1-naphthyl) ethylenediamine 2Hcl in a acid medium. The resulting color is pink. Any degree of pink color is considered positive.

pH: This test is based on double indicators (methyl red and bromothymol blue) which give a broad range of colors covering the entire urinary pH range. The color changes range from orange to greenish-yellow and green to blue.

Blood: This test is based on the pseudoperoxidase activity of hemoglobin which catalyses the reaction of O-Tolide and buffered organic peroxide and hydroperoxide. The resulting colors range from greenish-yellow to dark green.

Specific Gravity: This test is based on the pKa change of certain pre-treated polyelectrolytes in relation to ionic concentration. In the presence of an indicator, color changes from deep blue on urine with low ionic concentration to green and yellowish-green in urines of increasing ionic concentration.

Leucocytes: This test reveals the presence of granulocyte esterase on the leucocytes. These enzymes liberate indoxil ester and this indoxil reacts with a diazonium salt to produce a purple color.

KIT PRESENTATION

REF 524100-U - 100 determinations

1. Strips: 100 strips

2. Instructions for use

REF 524150-U - 150 determinations

1. Strips: 150 strips

2. Instructions for use

REF 524200-U - 200 determinations

1. Strips: 200 strips

2. Instructions for use

MATERIAL REQUIRED BUT NOT PROVIDED

STABILITY AND HANDLING

Store at 15 to 30°C. Do not use strips after the expiration date. Strips are sensitive to specific environmental factors such as moisture, heat and light; do not expose strips to these factors. Do not remove the desiccant from the bottle. All unused strips must be stored in the original bottle; therefore only strips to be tested should be removed from the tube. Darkening and discoloration of test area may indicate deterioration. If the deterioration is evident and the results are dubious, make sure the product is within the expiration date.

It is essential that strips are protected from humidity, heat and light.

SPECIMEN COLLECTION

Use a clean, dry, unused vessel to collect the urine. Fresh urine is recommended (up to 4 hours). Test urine as soon as possible after collection. If testing cannot be done within 4 hours, refrigerate the specimen immediately. Always label the specimen to reach room temperature before testing. Do not refrigerate the specimen, but do not centrifuge it.

PROCEDURE

1. Confirm if the products within the expiration date.

2. Remove the strip from the bottle and discard immediately.

3. Inspect carefully the strip, any signal of discoloration and darkening may indicate deterioration. In these cases, discard the strip.

4. Dip the test strip completely for 1 second in fresh, well-mixed and uncentrifuged urine specimen. Make sure all test areas are covered. Excessive urine may give wrong results. Remove the excessive urine by touching the plastic film on the rim of vessel. At this time, do not allow reagent areas to touch the rim of vessel. Excessive urine may be removed by gently blotting the lengthwise edge on a paper towel.

5. Compare the test results carefully with the color chart on the bottle label. Reading time from 30 to 60 seconds is critical for optimal results. While reading, keep the strip in a horizontal position to avoid possible interaction of chemicals by excessive urine. Changes in color that appear only along the edges of the test area or after more than 2 minutes have passed are of no diagnostic significance.

QUALITY CONTROL

The strips must be properly stored and handled before or during the testing. Reaction of reagent strips should be confirmed by testing known positive and negative specimens or multiple analytical controls containing normal and abnormal amounts of each of the analytes being tested.

RESULTS

The results are confirmed by direct comparison of test strip with the color charts printed on the bottle label. No calculations of laboratory instruments are necessary.

INTERPRETATION

Expected values: In this test strip, the normal urobilígen range is 0,1 to 1,0 mg/dL. If the results exceed the concentration of 2,0mg/dL, the patient and/or the urine specimen should be evaluated further.

Performance characteristics: This test can detect urobilígen in concentrations as low as 0,1 mg/dL; therefore, most normal urines may give a slight pink reaction.

Limitation of procedures: Higher concentration of formalin may give false negative result. This test is not reliable method for the detection of porphobilígen.

Chosen: Expected values: Amine quantity of glucose is excreted by the normal kidney. Concentration of 100mg/dL may be considered as abnormal if it is consistently found.

Performance characteristics: The test has a sensitivity of 50 to 100mg/dL and is specific for glucose. No substance excreted in urine other than glucose is known to give a positive result.

False negative results may be obtained with the presence of levodopa, ascorbic acid, glutathione and dipyrone. If the test color appears somewhat mottled at the higher glucose concentrations, match the darkest color to the color blocks.

Limitation of procedures: Reactivity of the test decreases as the specific gravity and/or pH of urine increases and may also vary with temperature. Ascorbic acid (more than 50mg/dL) and ketone bodies (more than 40mg/dL) may cause a false negative for a specimen containing a small amount of glucose (100mg/dL), however, the combinations of such ketones levels and low glucose levels are metabolically improbable in screening.

3. Ketones:

Expected values: Ketone bodies should not be detected in normal urine specimens with this reagent.

Performance characteristics: High specific gravity and low pH urines may give reactions up to and including trace level (5.0mg/dL).

Limitation of procedures: False positive results may occur with highly pigmented urine specimens or those containing large amounts of levodopa metabolites.

4. Bilirubin

Expected values: No bilirubin is detectable in urine of healthy people by even the most sensitive methods. Even trace amounts of bilirubin are sufficiently significant to require further investigation. Performance characteristics: The test has a sensitivity of 0.5mg/dl bilirubin. Limitation of procedures: Since bilirubin in specimens is sensitive to light, exposure of the urine specimen to light for a long period time may result in a false negative. Ascorbic acid concentration >25mg/dl may also cause false negative results. 5. Protein: Expected values: Normal urine specimens ordinarily contain low concentrations of protein (<20mg/dl); therefore, only persistent elevated levels of urine protein indicate kidney or urinary tract disease. The persistent results of trace level or over indicate significant proteinuria, and thus further clinical testing is needed to evaluate the significance of results. Limitation of procedures: The test is more sensitive to albumins than to gamma-globulins, Bence-Jones proteins, and mucoproteins; such proteins do not interfere with the reaction of albumin. Limitation of procedures: High buffered alkaline urine (pH 9) may give false positive results. The interpretation of results is also difficult in turbid urine specimens.

Expected values: No nitrite is detectable in urine specimens. Performance characteristics: The test is specific for nitrite and independent of urinary pH of any other substance normally excreted in urine. Comparison of the reacted area against a white background may aid in the detection of low levels of nitrite. Limitation of procedures: Ascobic acid concentration (>25mg/dl) can cause false negative on specimens with low nitrite concentration on urine (<0.03mg). False negative result should not be interpreted that patient is free of bacteriuria. Any degree of uniform pink color development should be considered positive; however, pink spots or pink edges should not be interpreted as a positive result. Negative results can occur when urinary tract infections are not caused by organisms with nitrate reductase, when the urine was not kept on the bladder for more than 4 hours to reduce nitrate to nitrite or when dietary nitrate is absent. 7. pH: Performance characteristics: Urine pH values generally range from 5 to 8 pH units. Performance characteristics: Values will be read within 1 unit over the pH value 5-9.

Limitation of procedures: Excessive urine on the test strip may wash the acid buffer from the neighboring protein reagent onto the pH area and change the pH to an acid pH although the urine being tested is originally neutral or alkaline. This is called the "run-over" phenomenon. 8. Blood/Hemoglobin: Expected values: No hemoglobin is detectable in urine (0.010mg/dl; 3 RBC/1). When hemoglobin appears in urine, it indicates some kind of kidney disease or some kind of urinary tract disorder. Blood is often, but not always, found in urine of menstruating females. Performance characteristics: This test is slightly more sensitive to free hemoglobin and myoglobin than the intact erythrocytes. The sensitivity may be reduced in urines with high specific gravity and ascorbic acid. The appearance of green spots on the reagent area indicates the presence of intact erythrocytes in the urine. Limitation of procedures: Ascobic acid or protein may reduce the reactivity of the blood test. Microbial peroxidase associated with urinary tract disorder can cause false positive results. Ascobic acid concentration (>40mg/dl) can cause false negative result on low level of blood. Urine from menstruating females may yield positive results. 9. Specific gravity: Expected values: The values range from 1.003 to 1.040. Performance characteristics: This test permits determination of urine specific gravity among 1,000-1,005-1,010-1,015-1,020-1,025-1,030. Highly buffered alkaline urines may cause low reading of result.

Limitation of procedures: Elevated specific gravity readings may be obtained in the presence of moderate quantities (100mg-700mg/dl) of protein. Expected values: Normally no leukocytes are detectable in urine. Performance characteristics: The test is generally capable of detecting 20-25WBC/ul as a trace. Microbial peroxidase associated with urinary tract disorder can cause false positive results. Ascobic acid concentration (>40mg/dl) can cause false negative result on low level of blood. Urine from menstruating females may yield positive results. High concentration of oxalic acid or trace of oxidizing agents may cause false negative results. 10. Leukocytes: 1. Urbilinogeno: 1-MeHyoxybenzenediazolium salt 2.5mg, Citrico acid 30mg 2. Glucose: Glucose oxidase 451 unit, Peroxidase 188 unit, Potassium iodide 100mg 3. Ketones: Sodium nitroprusside 24.00 mg, Magnesium sulfate 246.5mg 4. Bilirubin: 2,4-Dichlorobenzene-diazolium Na 0.3mg, Oxalic acid 30.0mg 5. Protein: Nitrobenzophenone blue 0.3mg, Citric acid 110.0mg, Sodium citrate 46.0mg 6. Nitrite: p-Arsanilic acid 5.0mg, N-(1-naphthyl) ethylenediamine 2Hd 6.0mg 7. pH: Methyl red 0.04mg, Bromothymol blue 0.5mg 8. Blood: Cumene hydroperoxide 7.00mg, O-Tolidine 3.00mg 9. Specific gravity: Bromothymol blue 1.2mg, Diethylenetriaminopentaacetic acid 12.0mg 10. Leukocytes: Phenylthiazole amino acid ester 1.0mg, Diazonium salt 0.7mg

SENSITIVITY AND SPECIFICITY Urbilinogeno: 0.1 mg/dl Glucose: 50-100mg/dl Ketones: 5mg/dl Bilirubin: 0.5mg/dl (bilirubin) Protein: 30mg/dl (albumin) Nitrite: 0.05mg/dl (nitrite ion) Blood: 10RBC/ul (0.03mg/dl hemoglobin intact erythrocyte) Leucocytes: 20-25WBC/ul (intact)

ANALYTICAL PRECAUTIONS In vitro diagnostic use only. This test must be used with urine specimens only. As all laboratory tests, final diagnostic or therapeutic decisions should not be based solely on one result or method. Effects of chemicals on the individual tests are not known in all cases. Therefore, it is recommended to repeat the test after removing such interfering agents such as drugs or vitamin supplements. Do not remove desiccant from tube. Do not touch the edge of a strip. Do not open the tube when ready for testing. The work area should be clean and free of detergents and other contaminants. Keep this product away from children. Use one test strip for each test. Follow the reading time printed on the label. Color changing on the edges of test area should be ignored; a careful removing of urine excess

eliminates the phenomenon. **PRECAUTIONS AND WARNINGS** 1. In vitro diagnostic use only. 2. Carefully read the package insert prior to testing. 3. Avoid exposing this product to direct sunlight. 4. Do not freeze the test strip since it can cause irreversible deterioration. 5. Store at 15-30°C. Do not use strips after the expiration date. 7. Dispose in accordance with local regulations. 8. Follow good manufacturing practices (GMP) on handling, storing and material disposal.

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1. Berg B, Helsing K, Jørgensen R, Kallner A. (1989). Guidelines for evaluation of reagent strips. Excerpted by analysis of urine albumin and glucose concentration using visually read reagent strips. International Federation of Clinical Chemistry, Scientific Division. **Scand J Clin Lab Invest.** 49(7):689-99. 2. Free AH, Free HM. (1986). Urinalysis: its proper role in the physician's office. **Clin Lab Med.** 1(2): 253-66. 3. Gyure WL. (1977). Comparison of several methods for semiquantitative determination of urinary protein. **Clin Chem.** 23(5):876-9. 4. Koss S, Porfi A, Wiedner A, Frank R, Vento S, Tractman H. (2006). Proteinuria and renal disease: prognocic value of urine dipstick testing for leukocytes. **Pediatr Nephrol.** 21(4):584-7. 5. Wilson LA. (2005). Urinalysis. **Nurs Stand.** 19(35):51-4.

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